

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PROPRA 10MG STPJ2A00
User ID :	pb97604	Batch NO. :	PJM2500610
Print Interval :	50000 Tablets	Running Time :	00 Hour 18 Min 21 Sec
Total Counter :	50114	Total Running Time :	2272 Hour 54 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	47	%	Manual ;	57	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		6400000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		35	Min	Supply Time ;	1.5	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	25.0 %	9.88	kN	No		
HEP :	20.0 %	9.49	kN	Yes	10	Rotations 20 Tablets
HCP :	2.0 %	8.07	kN	Yes	15	Times
Ref :		7.91	kN			
LCP :	2.0 %	7.75	kN	Yes	15	Times
LEP :	20.0 %	6.33	kN	No	15	Rotations 20 Tablets
LSP :	25.0 %	5.93	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	1.88	kN	Min :	1.86	kN	Punch NO.;	41
SD :	0.50	%	Max :	1.90	kN	Punch NO.;	11
Mean Main Pressure :	7.91	kN	Min :	6.66	kN	Punch NO.;	27
SD :	6.30	%	Max :	8.93	kN	Punch NO.;	10
Filling Depth :	Increase ;	13	Times	Decrease ;	23	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	195000	Tablets/hour
Feeder : Status ;	Auto		Speed ;	12	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	6.00	mm	Pre Thickness ;	3.00	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	4.44	mm	Main Thickness ;	1.44	mm
Filling Depth :	3.59	mm	Current Cam :	Cam B	
Main Air pressure :	0.813	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	5	Min			
Lower Punch Remain Time (S2):	34	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.32	kN		
	Max ;	0.35	kN	Max Punch ;	38
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.13	kN	Max Punch ;	38
Ejection Force :	Average ;	0.04	kN		
	Max ;	0.05	kN	Max Punch ;	11

Tablet Data (N/A)

Total Counter :	50114	Tablets			
A.W.C Counter :	40563	Tablets			
HEP :	18 (0.0%)	Tablets	LEP :	66 (0.2%)	Tablets
GOOD :	40479 (99.8%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PROPRA 10MG STPJ2A00
User ID :	pb97604	Batch NO. :	PJM2500610
Print Interval :	50000 Tablets	Running Time :	00 Hour 18 Min 21 Sec
Total Counter :	50168	Total Running Time :	2272 Hour 54 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	39	%	Manual ;	30	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		6400000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		35	Min	Supply Time ;	1.5	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	25.0 %	10.12	kN	No		
HEP :	20.0 %	9.71	kN	Yes	15	Rotations
HCP :	2.0 %	8.26	kN	Yes	20	Times
Ref :		8.10	kN			
LCP :	2.0 %	7.93	kN	Yes	20	Times
LEP :	20.0 %	6.48	kN	Yes	15	Rotations
LSP :	25.0 %	6.07	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	2.03	kN	Min :	1.99	kN	Punch NO.;	45
SD :	0.80	%	Max :	2.07	kN	Punch NO.;	1
Mean Main Pressure :	8.00	kN	Min :	6.93	kN	Punch NO.;	54
SD :	5.10	%	Max :	8.92	kN	Punch NO.;	10
Filling Depth :	Increase ;	2	Times	Decrease ;	4	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	195000	Tablets/hour
Feeder : Status ;	Auto		Speed ;	10	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	6.00	mm	Pre Thickness ;	3.00	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	4.64	mm	Main Thickness ;	1.64	mm
Filling Depth :	4.04	mm	Current Cam :	Cam B	
Main Air pressure :	0.813	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	5	Min			
Lower Punch Remain Time (S2):	34	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.32	kN		
	Max ;	0.35	kN	Max Punch ;	38
Lower Punch :	Average ;	0.07	kN		
	Max ;	0.13	kN	Max Punch ;	38
Ejection Force :	Average ;	0.04	kN		
	Max ;	0.05	kN	Max Punch ;	13

Tablet Data (N/A)

Total Counter :	50168	Tablets			
A.W.C Counter :	40563	Tablets			
HEP :	18 (0.0%)	Tablets	LEP :	66 (0.2%)	Tablets
GOOD :	40479 (99.8%)	Tablets			