

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PROPRA 10MG<br>STPJ2A00 |
| User ID :        | pb97604       | Batch NO. :          | PJM2500610              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 35 Min 26 Sec   |
| Total Counter :  | 300130        | Total Running Time : | 2274 Hour 11 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 48                | %       | Manual ;           | 57    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 6400000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 35                | Min     | Supply Time ;      | 1.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 40                | Min     | Supply Time ;      | 2.0   | Sec                  |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 25.0 % | 10.97             | kN      | No                 |       |                      |
| HEP :                              | 20.0 % | 10.53             | kN      | Yes                | 10    | Rotations 20 Tablets |
| HCP :                              | 2.0 %  | 8.95              | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 8.78              | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 8.60              | kN      | Yes                | 15    | Times                |
| LEP :                              | 20.0 % | 7.02              | kN      | No                 | 15    | Rotations 20 Tablets |
| LSP :                              | 25.0 % | 6.58              | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 1.83       | kN | Min : | 1.81       | kN | Punch NO.; | 21 |
| SD :                 | 0.50       | %  | Max : | 1.86       | kN | Punch NO.; | 1  |
| Mean Main Pressure : | 8.74       | kN | Min : | 7.26       | kN | Punch NO.; | 46 |
| SD :                 | 6.40       | %  | Max : | 9.85       | kN | Punch NO.; | 10 |
| Filling Depth :      | Increase ; | 16 | Times | Decrease ; | 39 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.00     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.40    | mm  | Main Thickness ; | 1.40     | mm           |
| Filling Depth :               | 3.56    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.816   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 5       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 34      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.30 | kN |             |    |
|                  | Max ;     | 0.35 | kN | Max Punch ; | 22 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 38 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 22 |

### Tablet Data (N/A)

|                 |                |         |       |            |         |
|-----------------|----------------|---------|-------|------------|---------|
| Total Counter : | 300130         | Tablets |       |            |         |
| A.W.C Counter : | 290046         | Tablets |       |            |         |
| HEP :           | 59 (0.0%)      | Tablets | LEP : | 263 (0.1%) | Tablets |
| GOOD :          | 289724 (99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PROPRA 10MG<br>STPJ2A00 |
| User ID :        | pb97604       | Batch NO. :          | PJM2500610              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 35 Min 26 Sec   |
| Total Counter :  | 300130        | Total Running Time : | 2274 Hour 11 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 46                | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 6400000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 35                | Min     | Supply Time ;      | 1.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 40                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 25.0 % | 10.95             | kN      | No                 |              |           |
| HEP :                              | 20.0 % | 10.52             | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %  | 8.94              | kN      | Yes 20             | Times        |           |
| Ref :                              |        | 8.76              | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 8.59              | kN      | Yes 20             | Times        |           |
| LEP :                              | 20.0 % | 7.01              | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 25.0 % | 6.57              | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.00       | kN | Min : | 1.97       | kN | Punch NO.; | 37 |
| SD :                 | 0.60       | %  | Max : | 2.04       | kN | Punch NO.; | 40 |
| Mean Main Pressure : | 8.86       | kN | Min : | 7.66       | kN | Punch NO.; | 54 |
| SD :                 | 4.90       | %  | Max : | 10.26      | kN | Punch NO.; | 10 |
| Filling Depth :      | Increase ; | 6  | Times | Decrease ; | 22 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12                | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.00              | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 4.60    | mm  | Main Thickness ; | 1.60              | mm           |
| Filling Depth :               | 4.03    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.816   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 5       | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 34      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.30 | kN |             |    |
|                  | Max ;     | 0.35 | kN | Max Punch ; | 22 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 38 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 22 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 300130            | Tablets |       |            |         |
| A.W.C Counter : | 290046            | Tablets |       |            |         |
| HEP :           | 59 (0.0%)         | Tablets | LEP : | 263 (0.1%) | Tablets |
| GOOD :          | 289724<br>(99.9%) | Tablets |       |            |         |