

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PROPRA 10MG<br>STPJ2A00 |
| User ID :        | rp98927       | Batch NO. :          | PJM2500610              |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 45 Min 09 Sec   |
| Total Counter :  | 1500810       | Total Running Time : | 2280 Hour 21 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 50        | %       | Manual ;           | 57           | RPM       |
| Filling Cam :                      |        | Cam B     |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |        | 4.5-10mm  |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |        | 6400000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 90        | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 190       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 60.0      | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 25.0 % | 11.25     | kN      | No                 |              |           |
| HEP :                              | 20.0 % | 10.80     | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 1.5 %  | 9.13      | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 9.00      | kN      |                    |              |           |
| LCP :                              | 1.5 %  | 8.86      | kN      | Yes 15             | Times        |           |
| LEP :                              | 20.0 % | 7.20      | kN      | No 15              | Rotations 20 | Tablets   |
| LSP :                              | 25.0 % | 6.75      | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 40.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None   |           |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |     |            |    |
|----------------------|------------|----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.61       | kN | Min : | 1.37       | kN  | Punch NO.; | 27 |
| SD :                 | 4.40       | %  | Max : | 1.73       | kN  | Punch NO.; | 6  |
| Mean Main Pressure : | 8.93       | kN | Min : | 5.80       | kN  | Punch NO.; | 27 |
| SD :                 | 10.40      | %  | Max : | 10.23      | kN  | Punch NO.; | 35 |
| Filling Depth :      | Increase ; | 58 | Times | Decrease ; | 114 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.55    | mm  | Pre Thickness ;  | 1.55     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.37    | mm  | Main Thickness ; | 1.37     | mm           |
| Filling Depth :               | 3.99    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.803   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 71      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 115     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.32 | kN |             |    |
|                  | Max ;     | 0.37 | kN | Max Punch ; | 13 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 24 |
| Ejection Force : | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 45 |

### Tablet Data (N/A)

|                 |                 |         |       |              |         |
|-----------------|-----------------|---------|-------|--------------|---------|
| Total Counter : | 1500810         | Tablets |       |              |         |
| A.W.C Counter : | 1484429         | Tablets |       |              |         |
| HEP :           | 656 (0.0%)      | Tablets | LEP : | 17286 (1.2%) | Tablets |
| GOOD :          | 1466487 (98.8%) | Tablets |       |              |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PROPRA 10MG<br>STPJ2A00 |
| User ID :        | rp98927       | Batch NO. :          | PJM2500610              |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 45 Min 09 Sec   |
| Total Counter :  | 1500810       | Total Running Time : | 2280 Hour 21 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 46                | %       | Manual ;           | 30    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 6400000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 90                | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 190               | Min     | Supply Time ;      | 1.0   | Sec                  |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 25.0 % | 11.27             | kN      | No                 |       |                      |
| HEP :                              | 20.0 % | 10.82             | kN      | Yes                | 15    | Rotations 20 Tablets |
| HCP :                              | 1.5 %  | 9.15              | kN      | Yes                | 20    | Times                |
| Ref :                              |        | 9.01              | kN      |                    |       |                      |
| LCP :                              | 1.5 %  | 8.88              | kN      | Yes                | 20    | Times                |
| LEP :                              | 20.0 % | 7.21              | kN      | Yes                | 15    | Rotations 20 Tablets |
| LSP :                              | 25.0 % | 6.76              | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |                   |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.85       | kN  | Min : | 1.76       | kN  | Punch NO.; | 27 |
| SD :                 | 2.30       | %   | Max : | 1.94       | kN  | Punch NO.; | 64 |
| Mean Main Pressure : | 9.12       | kN  | Min : | 7.66       | kN  | Punch NO.; | 55 |
| SD :                 | 5.50       | %   | Max : | 10.04      | kN  | Punch NO.; | 1  |
| Filling Depth :      | Increase ; | 154 | Times | Decrease ; | 174 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.55    | mm  | Pre Thickness ;  | 1.55     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.57    | mm  | Main Thickness ; | 1.57     | mm           |
| Filling Depth :               | 4.06    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.803   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 71      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 115     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.32 | kN |             |    |
|                  | Max ;     | 0.35 | kN | Max Punch ; | 52 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 24 |
| Ejection Force : | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 45 |

### Tablet Data (N/A)

|                 |                    |         |       |              |         |
|-----------------|--------------------|---------|-------|--------------|---------|
| Total Counter : | 1500810            | Tablets |       |              |         |
| A.W.C Counter : | 1484429            | Tablets |       |              |         |
| HEP :           | 656 (0.0%)         | Tablets | LEP : | 17286 (1.2%) | Tablets |
| GOOD :          | 1466487<br>(98.8%) | Tablets |       |              |         |