

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PROPRA 10MG<br>STPJ2A00 |
| User ID :        | sj97904       | Batch NO. :          | PJM2500610              |
| Print Interval : | 50000 Tablets | Running Time :       | 04 Hour 46 Min 56 Sec   |
| Total Counter :  | 5603160       | Total Running Time : | 2301 Hour 24 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 37                | %       | Manual ;           | 57           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 6400000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 90                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 190               | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 9.59              | kN      | No                 |              |           |
| HEP :                              | 22.0 %      | 9.00              | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 2.5 %       | 7.56              | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 7.38              | kN      |                    |              |           |
| LCP :                              | 2.5 %       | 7.19              | kN      | Yes 15             | Times        |           |
| LEP :                              | 28.0 %      | 5.31              | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 35.0 %      | 4.80              | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.80       | kN  | Min : | 1.62       | kN  | Punch NO.; | 26 |
| SD :                 | 3.60       | %   | Max : | 1.92       | kN  | Punch NO.; | 16 |
| Mean Main Pressure : | 7.43       | kN  | Min : | 5.86       | kN  | Punch NO.; | 27 |
| SD :                 | 8.00       | %   | Max : | 8.58       | kN  | Punch NO.; | 43 |
| Filling Depth :      | Increase ; | 155 | Times | Decrease ; | 223 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 9        | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.40    | mm  | Pre Thickness ;  | 1.40     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.37    | mm  | Main Thickness ; | 1.37     | mm           |
| Filling Depth :               | 3.55    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.749   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 87      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 161     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.32 | kN |             |    |
|                  | Max ;     | 0.37 | kN | Max Punch ; | 15 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 24 |
| Ejection Force : | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 62 |

### Tablet Data (N/A)

|                 |                 |         |       |              |         |
|-----------------|-----------------|---------|-------|--------------|---------|
| Total Counter : | 5603160         | Tablets |       |              |         |
| A.W.C Counter : | 5570486         | Tablets |       |              |         |
| HEP :           | 2217 (0.0%)     | Tablets | LEP : | 64461 (1.2%) | Tablets |
| GOOD :          | 5503808 (98.8%) | Tablets |       |              |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PROPRA 10MG<br>STPJ2A00 |
| User ID :        | sj97904       | Batch NO. :          | PJM2500610              |
| Print Interval : | 50000 Tablets | Running Time :       | 04 Hour 46 Min 56 Sec   |
| Total Counter :  | 5603160       | Total Running Time : | 2301 Hour 24 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 31                | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 6400000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 90                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 190               | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 25.0 % | 10.49             | kN      | No                 |              |           |
| HEP :                              | 20.0 % | 10.07             | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %  | 8.56              | kN      | Yes 20             | Times        |           |
| Ref :                              |        | 8.39              | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 8.22              | kN      | Yes 20             | Times        |           |
| LEP :                              | 20.0 % | 6.71              | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 25.0 % | 6.29              | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.80       | kN  | Min : | 1.74       | kN  | Punch NO.; | 22 |
| SD :                 | 1.50       | %   | Max : | 1.90       | kN  | Punch NO.; | 2  |
| Mean Main Pressure : | 8.44       | kN  | Min : | 7.49       | kN  | Punch NO.; | 22 |
| SD :                 | 4.60       | %   | Max : | 9.30       | kN  | Punch NO.; | 2  |
| Filling Depth :      | Increase ; | 322 | Times | Decrease ; | 350 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 8        | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.55    | mm  | Pre Thickness ;  | 1.55     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.57    | mm  | Main Thickness ; | 1.57     | mm           |
| Filling Depth :               | 4.04    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.749   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 87      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 161     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.32 | kN |             |    |
|                  | Max ;     | 0.37 | kN | Max Punch ; | 15 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 24 |
| Ejection Force : | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 62 |

### Tablet Data (N/A)

|                 |                 |         |       |              |         |
|-----------------|-----------------|---------|-------|--------------|---------|
| Total Counter : | 5603160         | Tablets |       |              |         |
| A.W.C Counter : | 5570551         | Tablets |       |              |         |
| HEP :           | 2217 (0.0%)     | Tablets | LEP : | 64461 (1.2%) | Tablets |
| GOOD :          | 5503873 (98.8%) | Tablets |       |              |         |