

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU BP400MG STIU2B00
User ID :	os190259	Batch NO. :	IUM2500940
Print Interval :	50000 Tablets	Running Time :	13 Hour 11 Min 44 Sec
Total Counter :	1991338	Total Running Time :	2317 Hour 54 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	100	%	Manual ;	57	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		2025000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	25.0 %	26.55	kN	No		
HEP :	20.0 %	25.49	kN	Yes 10	Rotations 20	Tablets
HCP :	2.0 %	21.67	kN	Yes 20	Times	
Ref :		21.24	kN			
LCP :	2.0 %	20.82	kN	Yes 20	Times	
LEP :	20.0 %	16.99	kN	Yes 15	Rotations 20	Tablets
LSP :	25.0 %	15.93	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	1	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.15	kN	Min :	1.94	kN	Punch NO.;	13
SD :	4.20	%	Max :	2.31	kN	Punch NO.;	3
Mean Main Pressure :	18.07	kN	Min :	16.26	kN	Punch NO.;	12
SD :	4.50	%	Max :	19.79	kN	Punch NO.;	34
Filling Depth :	Increase ;	67	Times	Decrease ;	42	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	22	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	22	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.80	mm	Pre Thickness ;	2.81	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	5.40	mm	Main Thickness ;	2.41	mm
Filling Depth :	10.38	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.843	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	5	Min			
Lower Punch Remain Time (S2):	17	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	On		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.05	kN		
	Max ;	0.17	kN	Max Punch ;	18
Lower Punch :	Average ;	0.02	kN		
	Max ;	0.07	kN	Max Punch ;	28
Ejection Force :	Average ;	0.05	kN		
	Max ;	0.09	kN	Max Punch ;	30

Tablet Data (N/A)

Total Counter :	1991338	Tablets			
A.W.C Counter :	1962771	Tablets			
HEP :	3453 (0.2%)	Tablets	LEP :	1166 (0.1%)	Tablets
GOOD :	1958152 (99.8%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU BP400MG STIU2B00
User ID :	os190259	Batch NO. :	IUM2500940
Print Interval :	50000 Tablets	Running Time :	13 Hour 11 Min 44 Sec
Total Counter :	1991338	Total Running Time :	2317 Hour 54 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	100	%	Manual ;	50	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		2025000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	25.0 %	26.80	kN	No		
HEP :	20.0 %	25.73	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	21.87	kN	Yes 20	Times	
Ref :		21.44	kN			
LCP :	2.0 %	21.01	kN	Yes 20	Times	
LEP :	20.0 %	17.15	kN	Yes 15	Rotations 20	Tablets
LSP :	25.0 %	16.08	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	2.64	kN	Min :	2.48	kN	Punch NO.;	13
SD :	3.10	%	Max :	2.87	kN	Punch NO.;	44
Mean Main Pressure :	20.79	kN	Min :	19.31	kN	Punch NO.;	30
SD :	3.40	%	Max :	22.84	kN	Punch NO.;	44
Filling Depth :	Increase ;	90	Times	Decrease ;	27	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	22	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	22	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.80	mm	Pre Thickness ;	2.81	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	5.40	mm	Main Thickness ;	2.41	mm
Filling Depth :	9.85	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.843	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	5	Min			
Lower Punch Remain Time (S2):	17	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	On		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.05	kN		
	Max ;	0.17	kN	Max Punch ;	18
Lower Punch :	Average ;	0.02	kN		
	Max ;	0.08	kN	Max Punch ;	28
Ejection Force :	Average ;	0.05	kN		
	Max ;	0.09	kN	Max Punch ;	30

Tablet Data (N/A)

Total Counter :	1991338	Tablets			
A.W.C Counter :	1962771	Tablets			
HEP :	3453 (0.2%)	Tablets	LEP :	1166 (0.1%)	Tablets
GOOD :	1958152 (99.8%)	Tablets			