

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Leveti1000mgSTLC2D00  |
| User ID :        | sj97904       | Batch NO. :          | LCM2503110            |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 03 Min 44 Sec |
| Total Counter :  | 100018        | Total Running Time : | 2386 Hour 43 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 168      | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 660000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 25.0 %      | 48.56    | kN      | No                 |              |           |
| HEP :                              | 20.0 %      | 46.61    | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 39.62    | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 38.84    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 38.07    | kN      | Yes 20             | Times        |           |
| LEP :                              | 20.0 %      | 31.08    | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 25.0 %      | 29.13    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 2.11  | kN    | Min :      | 2.06  | kN    | Punch NO.; | 19 |
|                      |            | 1.20  | %     | Max :      | 2.16  | kN    | Punch NO.; | 42 |
| Mean Main Pressure : |            | 39.46 | kN    | Min :      | 37.71 | kN    | Punch NO.; | 43 |
|                      | SD :       | 2.00  | %     | Max :      | 40.80 | kN    | Punch NO.; | 53 |
| Filling Depth :      | Increase ; | 49    | Times | Decrease ; | 60    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 25       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 9.00    | mm  | Pre Thickness ;  | 6.01     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.64    | mm  | Main Thickness ; | 3.65     | mm           |
| Filling Depth :               | 10.01   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.795   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 48      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 28      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | On       |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Lower Punch :    | Average ; | 0.15 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 32 |
| Ejection Force : | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 33 |

### Tablet Data (N/A)

|                 |                   |         |       |           |         |
|-----------------|-------------------|---------|-------|-----------|---------|
| Total Counter : | 100018            | Tablets |       |           |         |
| A.W.C Counter : | 98277             | Tablets |       |           |         |
| HEP :           | 0 (0%)            | Tablets | LEP : | 29 (0.0%) | Tablets |
| GOOD :          | 98248<br>(100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Leveti1000mgSTLC2D00  |
| User ID :        | sj97904       | Batch NO. :          | LCM2503110            |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 03 Min 44 Sec |
| Total Counter :  | 100018        | Total Running Time : | 2386 Hour 43 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 167      | %       | Manual ;           | 21           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 660000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 25.0 %      | 49.69    | kN      | No                 |              |           |
| HEP :                              | 20.0 %      | 47.70    | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 40.55    | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 39.75    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 38.96    | kN      | Yes 20             | Times        |           |
| LEP :                              | 20.0 %      | 31.80    | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 25.0 %      | 29.81    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 2.27  | kN    | Min :      | 2.22  | kN    | Punch NO.; | 2  |
|                      |            | 1.30  | %     | Max :      | 2.37  | kN    | Punch NO.; | 40 |
| Mean Main Pressure : |            | 39.95 | kN    | Min :      | 38.08 | kN    | Punch NO.; | 33 |
|                      | SD :       | 2.40  | %     | Max :      | 43.62 | kN    | Punch NO.; | 53 |
| Filling Depth :      | Increase ; | 34    | Times | Decrease ; | 61    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15     | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 25       | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 9.00   | mm  | Pre Thickness ;  | 6.01     | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.79   | mm  | Main Thickness ; | 3.80     | mm           |
| Filling Depth :               | 10.27  | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.795  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 48     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 28     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | On       |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Lower Punch :    | Average ; | 0.15 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 32 |
| Ejection Force : | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 33 |

### Tablet Data (N/A)

|                 |                   |         |       |           |         |
|-----------------|-------------------|---------|-------|-----------|---------|
| Total Counter : | 100018            | Tablets |       |           |         |
| A.W.C Counter : | 98277             | Tablets |       |           |         |
| HEP :           | 0 (0%)            | Tablets | LEP : | 29 (0.0%) | Tablets |
| GOOD :          | 98248<br>(100.0%) | Tablets |       |           |         |