

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	METFO 500mg STMR2A00
User ID :	vg96371	Batch NO. :	MRM2502750
Print Interval :	50000 Tablets	Running Time :	00 Hour 29 Min 46 Sec
Total Counter :	50054	Total Running Time :	2453 Hour 6 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	68	%	Manual ;	20	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		4500000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	25.0 %	12.80	kN	No		
HEP :	20.0 %	12.29	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	10.45	kN	Yes 15	Times	
Ref :		10.24	kN			
LCP :	2.0 %	10.04	kN	Yes 15	Times	
LEP :	20.0 %	8.19	kN	Yes 10	Rotations 15	Tablets
LSP :	25.0 %	7.68	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	1.25	kN	Min :	1.24	kN	Punch NO.;	38
SD :	0.60	%	Max :	1.28	kN	Punch NO.;	9
Mean Main Pressure :	10.18	kN	Min :	9.80	kN	Punch NO.;	43
SD :	2.20	%	Max :	11.03	kN	Punch NO.;	53
Filling Depth :	Increase ; 1	Times	Decrease ;	11	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	16	RPM	Capacity :	101760	Tablets/hour
Feeder : Status ;	Auto		Speed ;	11	RPM
Pre Upper Punch Position ;	2.00	mm			
Pre Lower Punch Position ;	6.70	mm	Pre Thickness ;	4.70	mm
Main Upper Punch Position ;	2.00	mm			
Main Lower Punch Position ;	4.95	mm	Main Thickness ;	2.95	mm
Filling Depth :	7.88	mm	Current Cam :	Cam C	
Main Air pressure :	0.822	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	11	Min			
Lower Punch Remain Time (S2):	16	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.07	kN		
	Max ;	0.14	kN	Max Punch ;	34
Lower Punch :	Average ;	0.11	kN		
	Max ;	0.16	kN	Max Punch ;	5
Ejection Force :	Average ;	1.97	kN		
	Max ;	1.99	kN	Max Punch ;	10

Tablet Data (N/A)

Total Counter :	50054	Tablets			
A.W.C Counter :	45569	Tablets			
HEP :	9 (0.0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	45560 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	METFO 500mg STMR2A00
User ID :	vg96371	Batch NO. :	MRM2502750
Print Interval :	50000 Tablets	Running Time :	00 Hour 29 Min 46 Sec
Total Counter :	50054	Total Running Time :	2453 Hour 6 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	66	%	Manual ;	21	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		4500000	Tablets			
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	25.0 %	12.18	kN	No		
HEP :	20.0 %	11.69	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	9.94	kN	Yes 15	Times	
Ref :		9.74	kN			
LCP :	2.0 %	9.55	kN	Yes 15	Times	
LEP :	20.0 %	7.80	kN	Yes 10	Rotations 15	Tablets
LSP :	25.0 %	7.31	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	1.36	kN	Min :	1.33	kN	Punch NO.;	8
SD :	0.70	%	Max :	1.38	kN	Punch NO.;	43
Mean Main Pressure :	9.82	kN	Min :	8.89	kN	Punch NO.;	28
SD :	3.40	%	Max :	10.71	kN	Punch NO.;	53
Filling Depth :	Increase ;	5	Times	Decrease ;	10	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	16	RPM	Capacity :	101760	Tablets/hour
Feeder : Status ;	Auto		Speed ;	11	RPM
Pre Upper Punch Position ;	2.00	mm			
Pre Lower Punch Position ;	6.80	mm	Pre Thickness ;	4.80	mm
Main Upper Punch Position ;	2.00	mm			
Main Lower Punch Position ;	5.10	mm	Main Thickness ;	3.10	mm
Filling Depth :	7.91	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.822	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	11	Min			
Lower Punch Remain Time (S2):	16	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.07	kN		
	Max ;	0.15	kN	Max Punch ;	34
Lower Punch :	Average ;	0.11	kN		
	Max ;	0.16	kN	Max Punch ;	5
Ejection Force :	Average ;	1.97	kN		
	Max ;	1.99	kN	Max Punch ;	10

Tablet Data (N/A)

Total Counter :	50054	Tablets			
A.W.C Counter :	45569	Tablets			
HEP :	9 (0.0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	45560 (100.0%)	Tablets			