

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU BP600MG STIU2C00
User ID :	rp98927	Batch NO. :	IUM2504460
Print Interval :	50000 Tablets	Running Time :	00 Hour 32 Min 03 Sec
Total Counter :	50032	Total Running Time :	2668 Hour 31 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	280	%	Manual ;	20	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		1350000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	5.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	30.83	kN	No		
HEP :	25.0 %	29.65	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	24.19	kN	Yes 15	Times	
Ref :		23.72	kN			
LCP :	2.0 %	23.24	kN	Yes 15	Times	
LEP :	25.0 %	17.79	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	16.60	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	4.02	kN	Min :	3.98	kN	Punch NO.;	15
SD :	0.60	%	Max :	4.08	kN	Punch NO.;	11
Mean Main Pressure :	23.92	kN	Min :	21.89	kN	Punch NO.;	37
SD :	4.40	%	Max :	26.57	kN	Punch NO.;	50
Filling Depth :	Increase ; 23	Times	Decrease ; 14	Times			

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	42	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	7.40	mm	Pre Thickness ;	4.40	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	7.29	mm	Main Thickness ;	4.29	mm
Filling Depth :	10.57	mm	Current Cam :	Cam C	
Main Air pressure :	0.824	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	46	Min			
Lower Punch Remain Time (S2):	6	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0
Lower Punch :	Average ;	0.07	kN		
	Max ;	0.08	kN	Max Punch ;	22
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	50032	Tablets			
A.W.C Counter :	43267	Tablets			
HEP :	47 (0.1%)	Tablets	LEP :	22 (0.1%)	Tablets
GOOD :	43198 (99.8%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU BP600MG STIU2C00
User ID :	rp98927	Batch NO. :	IUM2504460
Print Interval :	50000 Tablets	Running Time :	00 Hour 32 Min 03 Sec
Total Counter :	50032	Total Running Time :	2668 Hour 31 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	240	%	Manual ;	23	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		1350000	Tablets			
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	5.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	30.74	kN	No		
HEP :	25.0 %	29.55	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	24.12	kN	Yes 15	Times	
Ref :		23.64	kN			
LCP :	2.0 %	23.17	kN	Yes 15	Times	
LEP :	25.0 %	17.73	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	16.55	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	4.01	kN	Min :	3.97	kN	Punch NO.;	25
SD :	0.60	%	Max :	4.08	kN	Punch NO.;	53
Mean Main Pressure :	23.45	kN	Min :	20.92	kN	Punch NO.;	43
SD :	5.10	%	Max :	26.00	kN	Punch NO.;	27
Filling Depth :	Increase ;	10	Times	Decrease ;	13	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	36	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	7.40	mm	Pre Thickness ;	4.40	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	7.29	mm	Main Thickness ;	4.29	mm
Filling Depth :	9.55	mm	Current Cam :	Cam C	
Main Air pressure :	0.824	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	46	Min			
Lower Punch Remain Time (S2):	6	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0
Lower Punch :	Average ;	0.07	kN		
	Max ;	0.08	kN	Max Punch ;	22
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	50032	Tablets			
A.W.C Counter :	43267	Tablets			
HEP :	47 (0.1%)	Tablets	LEP :	22 (0.1%)	Tablets
GOOD :	43198 (99.8%)	Tablets			