

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Panto 40 mg STKQ1B00
User ID :	ms88737	Batch NO. :	KQM2509740
Print Interval :	50000 Tablets	Running Time :	01 Hour 40 Min 13 Sec
Total Counter :	324368	Total Running Time :	2867 Hour 27 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	60	%	Manual ;	23	RPM
Filling Cam :		Cam B		Cam Range :	Min :	2.50 mm
		4.5-10mm			Max :	8.00 mm
Target Quantity :		2750000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	34.18	kN	No		
HEP :	25.0 %	32.87	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	26.82	kN	Yes 15	Times	
Ref :		26.30	kN			
LCP :	2.0 %	25.77	kN	Yes 15	Times	
LEP :	25.0 %	19.72	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	18.41	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	SD :	3.76	kN	Min :	3.39	kN	Punch NO.;	20
		4.40	%	Max :	4.18	kN	Punch NO.;	65
Mean Main Pressure :		25.72	kN	Min :	20.21	kN	Punch NO.;	57
	SD :	6.20	%	Max :	29.67	kN	Punch NO.;	76
Filling Depth :	Increase ;	76	Times	Decrease ;	76	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	24	RPM	Capacity :	221760	Tablets/hour
Feeder : Status ;	Auto		Speed ;	14	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.20	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.90	mm	Main Thickness ;	0.90	mm
Filling Depth :	4.03	mm	Current Cam :	Cam B	
Main Air pressure :	0.730	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	11	Min			
Lower Punch Remain Time (S2):	6	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.23	kN		
	Max ;	0.36	kN	Max Punch ;	58
Lower Punch :	Average ;	0.19	kN		
	Max ;	0.54	kN	Max Punch ;	4
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	324368	Tablets			
A.W.C Counter :	288228	Tablets			
HEP :	145 (0.1%)	Tablets	LEP :	45 (0.0%)	Tablets
GOOD :	288038 (99.9%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Panto 40 mg STKQ1B00
User ID :	ms88737	Batch NO. :	KQM2509740
Print Interval :	50000 Tablets	Running Time :	01 Hour 40 Min 13 Sec
Total Counter :	324368	Total Running Time :	2867 Hour 27 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	60	%	Manual ;	27	RPM
Filling Cam :		Cam B		Cam Range :	Min :	2.50 mm
		4.5-10mm			Max :	8.00 mm
Target Quantity :		2750000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	35.95	kN	No		
HEP :	25.0 %	34.57	kN	Yes	10	Rotations 15 Tablets
HCP :	2.0 %	28.21	kN	Yes	15	Times
Ref :		27.65	kN			
LCP :	2.0 %	27.10	kN	Yes	15	Times
LEP :	25.0 %	20.74	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	19.36	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	SD :	3.98	kN	Min :	3.64	kN	Punch NO.;	52
		4.70	%	Max :	4.50	kN	Punch NO.;	2
Mean Main Pressure :		27.94	kN	Min :	22.99	kN	Punch NO.;	44
	SD :	5.40	%	Max :	30.72	kN	Punch NO.;	48
Filling Depth :	Increase ;	39	Times	Decrease ;	29	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	24	RPM	Capacity :	221760	Tablets/hour
Feeder : Status ;	Auto		Speed ;	14	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.20	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.84	mm	Main Thickness ;	0.84	mm
Filling Depth :	4.85	mm	Current Cam :	Cam B	
Main Air pressure :	0.730	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	11	Min			
Lower Punch Remain Time (S2):	6	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.23	kN		
	Max ;	0.33	kN	Max Punch ;	68
Lower Punch :	Average ;	0.19	kN		
	Max ;	0.54	kN	Max Punch ;	4
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	324368	Tablets			
A.W.C Counter :	288228	Tablets			
HEP :	145 (0.1%)	Tablets	LEP :	45 (0.0%)	Tablets
GOOD :	288038 (99.9%)	Tablets			