

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Panto 40 mg STKQ1B00  |
| User ID :        | rp98927       | Batch NO. :          | KQM2509740            |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 41 Min 58 Sec |
| Total Counter :  | 1675184       | Total Running Time : | 2874 Hour 29 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 60       | %       | Manual ;           | 23           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2750000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 50       | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 40.59    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 39.03    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 31.85    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 31.22    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 30.60    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 23.42    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 21.86    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 4.00           | kN    | Min :      | 3.70  | kN    | Punch NO.; | 28 |
| SD :                 | 3.80           | %     | Max :      | 4.52  | kN    | Punch NO.; | 10 |
| Mean Main Pressure : | 31.32          | kN    | Min :      | 26.13 | kN    | Punch NO.; | 44 |
| SD :                 | 4.60           | %     | Max :      | 35.37 | kN    | Punch NO.; | 62 |
| Filling Depth :      | Increase ; 201 | Times | Decrease ; | 218   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 184800   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.20    | mm  | Pre Thickness ;  | 1.20     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.84    | mm  | Main Thickness ; | 0.84     | mm           |
| Filling Depth :               | 4.12    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.788   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 17      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 12      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |   |
|------------------|-----------|------|----|-------------|---|
| Upper Punch :    | Average ; | 0.11 | kN |             |   |
|                  | Max ;     | 0.25 | kN | Max Punch ; | 1 |
| Lower Punch :    | Average ; | 0.16 | kN |             |   |
|                  | Max ;     | 0.48 | kN | Max Punch ; | 4 |
| Ejection Force : | Average ; | 0.00 | kN |             |   |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 1675184            | Tablets |       |            |         |
| A.W.C Counter : | 1634400            | Tablets |       |            |         |
| HEP :           | 273 (0.0%)         | Tablets | LEP : | 609 (0.0%) | Tablets |
| GOOD :          | 1633518<br>(99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Panto 40 mg STKQ1B00  |
| User ID :        | rp98927       | Batch NO. :          | KQM2509740            |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 41 Min 58 Sec |
| Total Counter :  | 1675184       | Total Running Time : | 2874 Hour 29 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 60       | %       | Manual ;           | 27           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2750000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 50       | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 34.76    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 33.43    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 27.28    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 26.74    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 26.21    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 20.06    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 18.72    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |               |       |            |       |       |            |    |
|----------------------|---------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 4.11          | kN    | Min :      | 3.77  | kN    | Punch NO.; | 5  |
| SD :                 | 4.00          | %     | Max :      | 4.67  | kN    | Punch NO.; | 72 |
| Mean Main Pressure : | 26.57         | kN    | Min :      | 22.71 | kN    | Punch NO.; | 44 |
| SD :                 | 4.70          | %     | Max :      | 31.15 | kN    | Punch NO.; | 34 |
| Filling Depth :      | Increase ; 95 | Times | Decrease ; | 89    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 184800   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.17    | mm  | Pre Thickness ;  | 1.17     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.90    | mm  | Main Thickness ; | 0.90     | mm           |
| Filling Depth :               | 4.88    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.788   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 17      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 12      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |   |
|------------------|-----------|------|----|-------------|---|
| Upper Punch :    | Average ; | 0.11 | kN |             |   |
|                  | Max ;     | 0.24 | kN | Max Punch ; | 1 |
| Lower Punch :    | Average ; | 0.16 | kN |             |   |
|                  | Max ;     | 0.48 | kN | Max Punch ; | 4 |
| Ejection Force : | Average ; | 0.00 | kN |             |   |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 1675184            | Tablets |       |            |         |
| A.W.C Counter : | 1634400            | Tablets |       |            |         |
| HEP :           | 273 (0.0%)         | Tablets | LEP : | 609 (0.0%) | Tablets |
| GOOD :          | 1633518<br>(99.9%) | Tablets |       |            |         |