

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Panto 40 mg STKQ1B00  |
| User ID :        | sj97904       | Batch NO. :          | KQM2509740            |
| Print Interval : | 50000 Tablets | Running Time :       | 10 Hour 36 Min 07 Sec |
| Total Counter :  | 2025422       | Total Running Time : | 2876 Hour 23 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 60       | %       | Manual ;           | 23           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2750000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 50       | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 40.59    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 39.03    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 31.85    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 31.22    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 30.60    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 23.42    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 21.86    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 4.09       | kN  | Min : | 3.62       | kN  | Punch NO.; | 5  |
| SD :                 | 4.50       | %   | Max : | 4.48       | kN  | Punch NO.; | 10 |
| Mean Main Pressure : | 31.85      | kN  | Min : | 25.67      | kN  | Punch NO.; | 44 |
| SD :                 | 4.80       | %   | Max : | 35.36      | kN  | Punch NO.; | 26 |
| Filling Depth :      | Increase ; | 271 | Times | Decrease ; | 280 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 184800            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12                | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 4.20    | mm  | Pre Thickness ;  | 1.20              | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 3.84    | mm  | Main Thickness ; | 0.84              | mm           |
| Filling Depth :               | 4.61    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.727   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 3       | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 18      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |   |
|------------------|-----------|------|----|-------------|---|
| Upper Punch :    | Average ; | 0.14 | kN |             |   |
|                  | Max ;     | 0.26 | kN | Max Punch ; | 1 |
| Lower Punch :    | Average ; | 0.15 | kN |             |   |
|                  | Max ;     | 0.47 | kN | Max Punch ; | 4 |
| Ejection Force : | Average ; | 0.00 | kN |             |   |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 2025422             | Tablets |       |            |         |
| A.W.C Counter : | 1983045             | Tablets |       |            |         |
| HEP :           | 274 (0.0%)          | Tablets | LEP : | 654 (0.0%) | Tablets |
| GOOD :          | 1982117<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Panto 40 mg STKQ1B00  |
| User ID :        | sj97904       | Batch NO. :          | KQM2509740            |
| Print Interval : | 50000 Tablets | Running Time :       | 10 Hour 36 Min 07 Sec |
| Total Counter :  | 2025422       | Total Running Time : | 2876 Hour 23 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |        |           |         |                    |       |                      |
|------------------------------------|--------|-----------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 60        | %       | Manual ;           | 27    | RPM                  |
| Filling Cam :                      |        | Cam B     |         | Cam Range :        | Min : | 2.50 mm              |
|                                    |        | 4.5-10mm  |         |                    | Max : | 8.00 mm              |
| Target Quantity :                  |        | 2750000   | Tablets |                    |       |                      |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |           |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 50        | Min     | Supply Time ;      | 0.2   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 40        | Min     | Supply Time ;      | 1.0   | Sec                  |
| Powder Supply Time :               |        | 60.0      | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 1         | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 2.00      | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 3.00      | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.90      | kN      |                    |       |                      |
|                                    |        |           |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 33.71     | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 32.42     | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 26.45     | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 25.93     | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 25.41     | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 19.45     | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 18.15     | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |           |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00     | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |           |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 3.91       | kN  | Min : | 3.57       | kN  | Punch NO.; | 5  |
| SD :                 | 3.90       | %   | Max : | 4.30       | kN  | Punch NO.; | 14 |
| Mean Main Pressure : | 25.74      | kN  | Min : | 21.72      | kN  | Punch NO.; | 44 |
| SD :                 | 5.10       | %   | Max : | 28.88      | kN  | Punch NO.; | 34 |
| Filling Depth :      | Increase ; | 109 | Times | Decrease ; | 101 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 184800   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.20    | mm  | Pre Thickness ;  | 1.20     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.90    | mm  | Main Thickness ; | 0.90     | mm           |
| Filling Depth :               | 4.94    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.727   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 3       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 18      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |   |
|------------------|-----------|------|----|-------------|---|
| Upper Punch :    | Average ; | 0.14 | kN |             |   |
|                  | Max ;     | 0.30 | kN | Max Punch ; | 1 |
| Lower Punch :    | Average ; | 0.15 | kN |             |   |
|                  | Max ;     | 0.47 | kN | Max Punch ; | 4 |
| Ejection Force : | Average ; | 0.00 | kN |             |   |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 2025422             | Tablets |       |            |         |
| A.W.C Counter : | 1983045             | Tablets |       |            |         |
| HEP :           | 274 (0.0%)          | Tablets | LEP : | 654 (0.0%) | Tablets |
| GOOD :          | 1982117<br>(100.0%) | Tablets |       |            |         |