

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Panto 40 mg STKQ1B00  |
| User ID :        | mr99799       | Batch NO. :          | KQM2510040            |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 42 Min 05 Sec |
| Total Counter :  | 748648        | Total Running Time : | 2911 Hour 39 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 42       | %       | Manual ;           | 23           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2750000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 40       | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 10.0         | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 38.62    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 37.13    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 30.30    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 29.70    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 29.11    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 22.28    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 20.79    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 4.06           | kN    | Min :      | 3.67  | kN    | Punch NO.; | 18 |
| SD :                 | 5.20           | %     | Max :      | 4.61  | kN    | Punch NO.; | 43 |
| Mean Main Pressure : | 29.26          | kN    | Min :      | 22.73 | kN    | Punch NO.; | 44 |
| SD :                 | 7.20           | %     | Max :      | 33.84 | kN    | Punch NO.; | 5  |
| Filling Depth :      | Increase ; 117 | Times | Decrease ; | 102   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 203280   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 9        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.20    | mm  | Pre Thickness ;  | 1.21     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.81    | mm  | Main Thickness ; | 0.82     | mm           |
| Filling Depth :               | 3.75    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.835   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 7       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 18      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.26 | kN |             |    |
|                  | Max ;     | 0.39 | kN | Max Punch ; | 62 |
| Lower Punch :    | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.53 | kN | Max Punch ; | 4  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |             |         |
|-----------------|-------------------|---------|-------|-------------|---------|
| Total Counter : | 748648            | Tablets |       |             |         |
| A.W.C Counter : | 739978            | Tablets |       |             |         |
| HEP :           | 162 (0.0%)        | Tablets | LEP : | 1107 (0.1%) | Tablets |
| GOOD :          | 738709<br>(99.8%) | Tablets |       |             |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Panto 40 mg STKQ1B00  |
| User ID :        | mr99799       | Batch NO. :          | KQM2510040            |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 42 Min 05 Sec |
| Total Counter :  | 748648        | Total Running Time : | 2911 Hour 39 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 43       | %       | Manual ;           | 27           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2750000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 40       | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 10.0         | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 35.97    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 34.59    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 28.23    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 27.67    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 27.12    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 20.75    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 19.37    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |               |       |               |       |    |            |    |
|----------------------|---------------|-------|---------------|-------|----|------------|----|
| Mean Pre Pressure :  | 4.08          | kN    | Min :         | 3.66  | kN | Punch NO.; | 5  |
| SD :                 | 4.10          | %     | Max :         | 4.50  | kN | Punch NO.; | 72 |
| Mean Main Pressure : | 28.05         | kN    | Min :         | 23.88 | kN | Punch NO.; | 71 |
| SD :                 | 5.70          | %     | Max :         | 33.39 | kN | Punch NO.; | 40 |
| Filling Depth :      | Increase ; 74 | Times | Decrease ; 71 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 203280   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 9        | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.20    | mm  | Pre Thickness ;  | 1.20     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.86    | mm  | Main Thickness ; | 0.87     | mm           |
| Filling Depth :               | 3.72    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.835   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 7       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 18      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.26 | kN |             |    |
|                  | Max ;     | 0.36 | kN | Max Punch ; | 26 |
| Lower Punch :    | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.53 | kN | Max Punch ; | 4  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |             |         |
|-----------------|-------------------|---------|-------|-------------|---------|
| Total Counter : | 748648            | Tablets |       |             |         |
| A.W.C Counter : | 739978            | Tablets |       |             |         |
| HEP :           | 162 (0.0%)        | Tablets | LEP : | 1107 (0.1%) | Tablets |
| GOOD :          | 738709<br>(99.8%) | Tablets |       |             |         |