

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	sj97904	Batch NO. :	KQM2510240
Print Interval :	50000 Tablets	Running Time :	14 Hour 24 Min 22 Sec
Total Counter :	2676644	Total Running Time :	2960 Hour 37 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	70	%	Manual ;	23	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	40.40	kN	No		
HEP :	25.0 %	38.85	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	31.70	kN	Yes 15	Times	
Ref :		31.08	kN			
LCP :	2.0 %	30.46	kN	Yes 15	Times	
LEP :	25.0 %	23.31	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	21.75	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	SD :	3.76	kN	Min :	3.33	kN	Punch NO.;	65
		4.80	%	Max :	4.11	kN	Punch NO.;	33
Mean Main Pressure :	SD :	30.46	kN	Min :	27.29	kN	Punch NO.;	66
		4.70	%	Max :	33.40	kN	Punch NO.;	51
Filling Depth :	Increase ;	386	Times	Decrease ;	380	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	184800	Tablets/hour
Feeder : Status ;	Auto		Speed ;	14	RPM
Pre Upper Punch Position ;	2.69	mm			
Pre Lower Punch Position ;	4.00	mm	Pre Thickness ;	1.31	mm
Main Upper Punch Position ;	2.69	mm			
Main Lower Punch Position ;	3.60	mm	Main Thickness ;	0.91	mm
Filling Depth :	4.51	mm	Current Cam :	Cam B	
Main Air pressure :	0.724	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	19	Min			
Lower Punch Remain Time (S2):	14	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.11	kN		
	Max ;	0.18	kN	Max Punch ;	32
Lower Punch :	Average ;	0.38	kN		
	Max ;	0.82	kN	Max Punch ;	9
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	2676644	Tablets			
A.W.C Counter :	2636585	Tablets			
HEP :	443 (0.0%)	Tablets	LEP :	2012 (0.1%)	Tablets
GOOD :	2634130 (99.9%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	sj97904	Batch NO. :	KQM2510240
Print Interval :	50000 Tablets	Running Time :	14 Hour 24 Min 22 Sec
Total Counter :	2676644	Total Running Time :	2960 Hour 37 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	59	%	Manual ;	24	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	38.38	kN	No		
HEP :	25.0 %	36.91	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	30.12	kN	Yes 15	Times	
Ref :		29.52	kN			
LCP :	2.0 %	28.93	kN	Yes 15	Times	
LEP :	25.0 %	22.14	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	20.67	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	3.82	kN	Min :	3.39	kN	Punch NO.;	77
SD :	5.30	%	Max :	4.31	kN	Punch NO.;	38
Mean Main Pressure :	29.87	kN	Min :	26.47	kN	Punch NO.;	77
SD :	4.60	%	Max :	32.82	kN	Punch NO.;	48
Filling Depth :	Increase ;	438	Times	Decrease ;	438	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	184800	Tablets/hour
Feeder : Status ;	Auto		Speed ;	12	RPM
Pre Upper Punch Position ;	2.70	mm			
Pre Lower Punch Position ;	4.00	mm	Pre Thickness ;	1.30	mm
Main Upper Punch Position ;	2.75	mm			
Main Lower Punch Position ;	3.60	mm	Main Thickness ;	0.85	mm
Filling Depth :	4.10	mm	Current Cam :	Cam B	
Main Air pressure :	0.724	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	19	Min			
Lower Punch Remain Time (S2):	14	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.10	kN		
	Max ;	0.17	kN	Max Punch ;	32
Lower Punch :	Average ;	0.38	kN		
	Max ;	0.82	kN	Max Punch ;	9
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	2676644	Tablets			
A.W.C Counter :	2636585	Tablets			
HEP :	443 (0.0%)	Tablets	LEP :	2012 (0.1%)	Tablets
GOOD :	2634130 (99.9%)	Tablets			