

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	sj97904	Batch NO. :	KQM2510240
Print Interval :	50000 Tablets	Running Time :	14 Hour 43 Min 42 Sec
Total Counter :	2735996	Total Running Time :	2960 Hour 56 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	62	%	Manual ;	23	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	40.04	kN	No		
HEP :	25.0 %	38.50	kN	Yes	10	Rotations 15 Tablets
HCP :	2.0 %	31.42	kN	Yes	15	Times
Ref :		30.80	kN			
LCP :	2.0 %	30.18	kN	Yes	15	Times
LEP :	25.0 %	23.10	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	21.56	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	3.96	kN	Min :	3.46	kN	Punch NO.;	19
SD :	9.20	%	Max :	5.34	kN	Punch NO.;	72
Mean Main Pressure :	31.52	kN	Min :	27.85	kN	Punch NO.;	19
SD :	8.30	%	Max :	40.21	kN	Punch NO.;	72
Filling Depth :	Increase ;	424	Times	Decrease ;	398	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	12	RPM
Pre Upper Punch Position ;	2.69	mm			
Pre Lower Punch Position ;	4.00	mm	Pre Thickness ;	1.31	mm
Main Upper Punch Position ;	2.69	mm			
Main Lower Punch Position ;	3.60	mm	Main Thickness ;	0.91	mm
Filling Depth :	4.45	mm	Current Cam :	Cam B 4.5-10mm	
Main Air pressure :	0.747	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	1	Min			
Lower Punch Remain Time (S2):	16	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	On		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.11	kN		
	Max ;	0.17	kN	Max Punch ;	29
Lower Punch :	Average ;	0.40	kN		
	Max ;	0.79	kN	Max Punch ;	49
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	2735996	Tablets			
A.W.C Counter :	2695056	Tablets			
HEP :	461 (0.0%)	Tablets	LEP :	2384 (0.1%)	Tablets
GOOD :	2692211 (99.9%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	sj97904	Batch NO. :	KQM2510240
Print Interval :	50000 Tablets	Running Time :	14 Hour 43 Min 42 Sec
Total Counter :	2735996	Total Running Time :	2960 Hour 56 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	59	%	Manual ;	24	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	38.38	kN	No		
HEP :	25.0 %	36.91	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	30.12	kN	Yes 15	Times	
Ref :		29.52	kN			
LCP :	2.0 %	28.93	kN	Yes 15	Times	
LEP :	25.0 %	22.14	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	20.67	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	4.13	kN	Min :	3.59	kN	Punch NO.;	63
SD :	12.00	%	Max :	5.99	kN	Punch NO.;	32
Mean Main Pressure :	31.28	kN	Min :	27.69	kN	Punch NO.;	63
SD :	8.60	%	Max :	40.82	kN	Punch NO.;	32
Filling Depth :	Increase ;	465	Times	Decrease ;	451	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	12	RPM
Pre Upper Punch Position ;	2.70	mm			
Pre Lower Punch Position ;	4.00	mm	Pre Thickness ;	1.30	mm
Main Upper Punch Position ;	2.75	mm			
Main Lower Punch Position ;	3.60	mm	Main Thickness ;	0.85	mm
Filling Depth :	4.52	mm	Current Cam :	Cam B 4.5-10mm	
Main Air pressure :	0.747	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	1	Min			
Lower Punch Remain Time (S2):	16	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	On		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.12	kN		
	Max ;	0.37	kN	Max Punch ;	58
Lower Punch :	Average ;	0.40	kN		
	Max ;	0.79	kN	Max Punch ;	49
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	2735996	Tablets			
A.W.C Counter :	2695056	Tablets			
HEP :	461 (0.0%)	Tablets	LEP :	2384 (0.1%)	Tablets
GOOD :	2692211 (99.9%)	Tablets			