

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40<br>MG,STKQ1B00 |
| User ID :        | pb97604       | Batch NO. :          | KQM2510140              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 54 Min 58 Sec   |
| Total Counter :  | 350244        | Total Running Time : | 2962 Hour 51 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 62                | %       | Manual ;           | 23           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 2750000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 20                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 10.0         | Sec       |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 39.77             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 38.24             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 31.21             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 30.60             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 29.98             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 22.95             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 21.42             | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 3.77  | kN    | Min :      | 3.37  | kN    | Punch NO.; | 41 |
|                      |            | 3.50  | %     | Max :      | 4.08  | kN    | Punch NO.; | 74 |
| Mean Main Pressure : | SD :       | 30.78 | kN    | Min :      | 27.45 | kN    | Punch NO.; | 37 |
|                      |            | 3.10  | %     | Max :      | 32.72 | kN    | Punch NO.; | 74 |
| Filling Depth :      | Increase ; | 27    | Times | Decrease ; | 22    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12                | RPM          |
| Pre Upper Punch Position ;    | 2.69    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 4.00    | mm  | Pre Thickness ;  | 1.31              | mm           |
| Main Upper Punch Position ;   | 2.69    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 3.60    | mm  | Main Thickness ; | 0.91              | mm           |
| Filling Depth :               | 3.94    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.763   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 1       | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 16      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 12 |
| Lower Punch :    | Average ; | 0.19 | kN |             |    |
|                  | Max ;     | 0.48 | kN | Max Punch ; | 9  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 350244            | Tablets |       |            |         |
| A.W.C Counter : | 346490            | Tablets |       |            |         |
| HEP :           | 0 (0%)            | Tablets | LEP : | 308 (0.1%) | Tablets |
| GOOD :          | 346182<br>(99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40<br>MG,STKQ1B00 |
| User ID :        | pb97604       | Batch NO. :          | KQM2510140              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 54 Min 58 Sec   |
| Total Counter :  | 350244        | Total Running Time : | 2962 Hour 51 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 59                | %       | Manual ;           | 24           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 2750000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 20                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 10.0         | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.90              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 37.58             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 36.13             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 29.48             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 28.91             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 28.33             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 21.68             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 20.23             | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 1      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 3.65       | kN | Min : | 3.32       | kN | Punch NO.; | 5  |
| SD :                 | 4.70       | %  | Max : | 4.12       | kN | Punch NO.; | 36 |
| Mean Main Pressure : | 29.17      | kN | Min : | 26.54      | kN | Punch NO.; | 15 |
| SD :                 | 4.00       | %  | Max : | 31.77      | kN | Punch NO.; | 36 |
| Filling Depth :      | Increase ; | 34 | Times | Decrease ; | 30 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12                | RPM          |
| Pre Upper Punch Position ;    | 2.70    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 4.00    | mm  | Pre Thickness ;  | 1.30              | mm           |
| Main Upper Punch Position ;   | 2.75    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 3.60    | mm  | Main Thickness ; | 0.85              | mm           |
| Filling Depth :               | 3.36    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.763   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 1       | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 16      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 12 |
| Lower Punch :    | Average ; | 0.20 | kN |             |    |
|                  | Max ;     | 0.46 | kN | Max Punch ; | 9  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 350244            | Tablets |       |            |         |
| A.W.C Counter : | 346490            | Tablets |       |            |         |
| HEP :           | 0 (0%)            | Tablets | LEP : | 308 (0.1%) | Tablets |
| GOOD :          | 346182<br>(99.9%) | Tablets |       |            |         |