

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	sj97904	Batch NO. :	KQM2510340
Print Interval :	50000 Tablets	Running Time :	09 Hour 48 Min 49 Sec
Total Counter :	1805756	Total Running Time :	2985 Hour 28 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	54	%	Manual ;	21	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	38.13	kN	No		
HEP :	25.0 %	36.66	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	29.92	kN	Yes 15	Times	
Ref :		29.33	kN			
LCP :	2.0 %	28.74	kN	Yes 15	Times	
LEP :	25.0 %	22.00	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	20.53	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	3.82	kN	Min :	3.50	kN	Punch NO.;	71
SD :	4.50	%	Max :	4.34	kN	Punch NO.;	48
Mean Main Pressure :	29.07	kN	Min :	26.44	kN	Punch NO.;	71
SD :	4.10	%	Max :	32.45	kN	Punch NO.;	48
Filling Depth :	Increase ;	252	Times	Decrease ;	192	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	184800	Tablets/hour
Feeder : Status ;	Auto		Speed ;	11	RPM
Pre Upper Punch Position ;	2.69	mm			
Pre Lower Punch Position ;	4.00	mm	Pre Thickness ;	1.31	mm
Main Upper Punch Position ;	2.69	mm			
Main Lower Punch Position ;	3.68	mm	Main Thickness ;	0.99	mm
Filling Depth :	3.71	mm	Current Cam :	Cam B	
Main Air pressure :	0.765	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	45	Min			
Lower Punch Remain Time (S2):	13	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.17	kN		
	Max ;	0.22	kN	Max Punch ;	12
Lower Punch :	Average ;	0.17	kN		
	Max ;	0.40	kN	Max Punch ;	9
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	1805756	Tablets			
A.W.C Counter :	1795478	Tablets			
HEP :	335 (0.0%)	Tablets	LEP :	469 (0.0%)	Tablets
GOOD :	1794674 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	sj97904	Batch NO. :	KQM2510340
Print Interval :	50000 Tablets	Running Time :	09 Hour 48 Min 49 Sec
Total Counter :	1805756	Total Running Time :	2985 Hour 28 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	52	%	Manual ;	24	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	36.58	kN	No		
HEP :	25.0 %	35.17	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	28.70	kN	Yes 15	Times	
Ref :		28.14	kN			
LCP :	2.0 %	27.58	kN	Yes 15	Times	
LEP :	25.0 %	21.10	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	19.70	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	3.73	kN	Min :	3.30	kN	Punch NO.;	9
SD :	6.10	%	Max :	4.33	kN	Punch NO.;	43
Mean Main Pressure :	27.68	kN	Min :	23.92	kN	Punch NO.;	9
SD :	5.20	%	Max :	30.93	kN	Punch NO.;	43
Filling Depth :	Increase ;	515	Times	Decrease ;	512	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	184800	Tablets/hour
Feeder : Status ;	Auto		Speed ;	10	RPM
Pre Upper Punch Position ;	2.70	mm			
Pre Lower Punch Position ;	4.00	mm	Pre Thickness ;	1.30	mm
Main Upper Punch Position ;	2.74	mm			
Main Lower Punch Position ;	3.68	mm	Main Thickness ;	0.94	mm
Filling Depth :	3.54	mm	Current Cam :	Cam B	
Main Air pressure :	0.765	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	45	Min			
Lower Punch Remain Time (S2):	13	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.17	kN		
	Max ;	0.22	kN	Max Punch ;	12
Lower Punch :	Average ;	0.17	kN		
	Max ;	0.41	kN	Max Punch ;	9
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	1805756	Tablets			
A.W.C Counter :	1795478	Tablets			
HEP :	335 (0.0%)	Tablets	LEP :	469 (0.0%)	Tablets
GOOD :	1794674 (100.0%)	Tablets			