

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	sj97904	Batch NO. :	KQM2510440
Print Interval :	50000 Tablets	Running Time :	03 Hour 32 Min 50 Sec
Total Counter :	650306	Total Running Time :	2994 Hour 0 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	68	%	Manual ;	21	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	46.48	kN	No		
HEP :	25.0 %	44.69	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	36.47	kN	Yes 15	Times	
Ref :		35.76	kN			
LCP :	2.0 %	35.04	kN	Yes 15	Times	
LEP :	25.0 %	26.82	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	25.03	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	3.85	kN	Min :	3.30	kN	Punch NO.;	17
SD :	7.00	%	Max :	4.53	kN	Punch NO.;	36
Mean Main Pressure :	35.14	kN	Min :	31.46	kN	Punch NO.;	16
SD :	5.20	%	Max :	39.65	kN	Punch NO.;	36
Filling Depth :	Increase ;	272	Times	Decrease ;	256	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	14	RPM
Pre Upper Punch Position ;	2.70	mm			
Pre Lower Punch Position ;	4.00	mm	Pre Thickness ;	1.30	mm
Main Upper Punch Position ;	2.70	mm			
Main Lower Punch Position ;	3.52	mm	Main Thickness ;	0.82	mm
Filling Depth :	4.42	mm	Current Cam :	Cam B 4.5-10mm	
Main Air pressure :	0.754	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	3	Min			
Lower Punch Remain Time (S2):	15	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.21	kN		
	Max ;	0.29	kN	Max Punch ;	61
Lower Punch :	Average ;	0.46	kN		
	Max ;	1.07	kN	Max Punch ;	46
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	650306	Tablets			
A.W.C Counter :	635370	Tablets			
HEP :	54 (0.0%)	Tablets	LEP :	231 (0.0%)	Tablets
GOOD :	635085 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	sj97904	Batch NO. :	KQM2510440
Print Interval :	50000 Tablets	Running Time :	03 Hour 32 Min 50 Sec
Total Counter :	650306	Total Running Time :	2994 Hour 0 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	53	%	Manual ;	24	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	39.95	kN	No		
HEP :	25.0 %	38.41	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	31.35	kN	Yes 15	Times	
Ref :		30.73	kN			
LCP :	2.0 %	30.12	kN	Yes 15	Times	
LEP :	25.0 %	23.05	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	21.51	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	3.72	kN	Min :	3.41	kN	Punch NO.;	17
SD :	3.90	%	Max :	4.03	kN	Punch NO.;	32
Mean Main Pressure :	31.14	kN	Min :	28.44	kN	Punch NO.;	6
SD :	3.30	%	Max :	33.36	kN	Punch NO.;	32
Filling Depth :	Increase ;	56	Times	Decrease ;	61	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	11	RPM
Pre Upper Punch Position ;	2.70	mm			
Pre Lower Punch Position ;	4.00	mm	Pre Thickness ;	1.30	mm
Main Upper Punch Position ;	2.70	mm			
Main Lower Punch Position ;	3.50	mm	Main Thickness ;	0.80	mm
Filling Depth :	3.47	mm	Current Cam :	Cam B 4.5-10mm	
Main Air pressure :	0.754	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	3	Min			
Lower Punch Remain Time (S2):	15	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.21	kN		
	Max ;	0.29	kN	Max Punch ;	61
Lower Punch :	Average ;	0.46	kN		
	Max ;	1.07	kN	Max Punch ;	46
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	650306	Tablets			
A.W.C Counter :	635370	Tablets			
HEP :	54 (0.0%)	Tablets	LEP :	231 (0.0%)	Tablets
GOOD :	635085 (100.0%)	Tablets			