

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	sp98936	Batch NO. :	KQM2510440
Print Interval :	50000 Tablets	Running Time :	07 Hour 59 Min 34 Sec
Total Counter :	1500728	Total Running Time :	2998 Hour 27 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	62	%	Manual ;	45	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	37.07	kN	No		
HEP :	25.0 %	35.64	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	29.08	kN	Yes 15	Times	
Ref :		28.51	kN			
LCP :	2.0 %	27.94	kN	Yes 15	Times	
LEP :	25.0 %	21.38	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	19.96	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	3.90	kN	Min :	3.47	kN	Punch NO.;	52
SD :	6.10	%	Max :	4.59	kN	Punch NO.;	77
Mean Main Pressure :	28.38	kN	Min :	25.48	kN	Punch NO.;	65
SD :	5.20	%	Max :	31.75	kN	Punch NO.;	77
Filling Depth :	Increase ;	486	Times	Decrease ;	465	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	12	RPM
Pre Upper Punch Position ;	2.70	mm			
Pre Lower Punch Position ;	3.90	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	2.70	mm			
Main Lower Punch Position ;	3.55	mm	Main Thickness ;	0.85	mm
Filling Depth :	4.60	mm	Current Cam :	Cam B	
Main Air pressure :	0.792	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	48	Min			
Lower Punch Remain Time (S2):	20	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	On		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.14	kN		
	Max ;	0.21	kN	Max Punch ;	21
Lower Punch :	Average ;	0.16	kN		
	Max ;	0.47	kN	Max Punch ;	72
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	1500728	Tablets			
A.W.C Counter :	1461351	Tablets			
HEP :	286 (0.0%)	Tablets	LEP :	787 (0.1%)	Tablets
GOOD :	1460278 (99.9%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	sp98936	Batch NO. :	KQM2510440
Print Interval :	50000 Tablets	Running Time :	07 Hour 59 Min 34 Sec
Total Counter :	1500728	Total Running Time :	2998 Hour 27 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	64	%	Manual ;	45	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	33.82	kN	No		
HEP :	25.0 %	32.52	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	26.54	kN	Yes 15	Times	
Ref :		26.02	kN			
LCP :	2.0 %	25.50	kN	Yes 15	Times	
LEP :	25.0 %	19.51	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	18.21	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	4.18	kN	Min :	3.50	kN	Punch NO.;	18
SD :	9.80	%	Max :	5.53	kN	Punch NO.;	76
Mean Main Pressure :	28.04	kN	Min :	23.93	kN	Punch NO.;	18
SD :	6.60	%	Max :	34.34	kN	Punch NO.;	76
Filling Depth :	Increase ;	179	Times	Decrease ;	152	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	13	RPM
Pre Upper Punch Position ;	2.70	mm			
Pre Lower Punch Position ;	3.90	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	2.70	mm			
Main Lower Punch Position ;	3.55	mm	Main Thickness ;	0.85	mm
Filling Depth :	4.20	mm	Current Cam :	Cam B	
Main Air pressure :	0.792	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	48	Min			
Lower Punch Remain Time (S2):	20	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	On		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.16	kN		
	Max ;	0.30	kN	Max Punch ;	74
Lower Punch :	Average ;	0.16	kN		
	Max ;	0.47	kN	Max Punch ;	72
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	1500728	Tablets			
A.W.C Counter :	1461351	Tablets			
HEP :	286 (0.0%)	Tablets	LEP :	787 (0.1%)	Tablets
GOOD :	1460278 (99.9%)	Tablets			