

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	pb97604	Batch NO. :	KQM2510440
Print Interval :	50000 Tablets	Running Time :	08 Hour 56 Min 18 Sec
Total Counter :	1700856	Total Running Time :	2999 Hour 24 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	74	%	Manual ;	45	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	40.91	kN	No		
HEP :	25.0 %	39.34	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	32.10	kN	Yes 15	Times	
Ref :		31.47	kN			
LCP :	2.0 %	30.84	kN	Yes 15	Times	
LEP :	25.0 %	23.60	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	22.03	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	4.00	kN	Min :	3.53	kN	Punch NO.;	62
SD :	5.80	%	Max :	4.53	kN	Punch NO.;	11
Mean Main Pressure :	30.75	kN	Min :	27.54	kN	Punch NO.;	41
SD :	4.90	%	Max :	34.30	kN	Punch NO.;	11
Filling Depth :	Increase ; 552	Times	Decrease ;	512	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	24	RPM	Capacity :	221760	Tablets/hour
Feeder : Status ;	Auto		Speed ;	18	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.21	mm	Pre Thickness ;	1.21	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.84	mm	Main Thickness ;	0.84	mm
Filling Depth :	4.79	mm	Current Cam :	Cam B	
Main Air pressure :	0.642	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	45	Min			
Lower Punch Remain Time (S2):	7	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.62	kN		
	Max ;	0.85	kN	Max Punch ;	8
Lower Punch :	Average ;	0.44	kN		
	Max ;	1.02	kN	Max Punch ;	49
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	1700856	Tablets			
A.W.C Counter :	1657793	Tablets			
HEP :	334 (0.0%)	Tablets	LEP :	1028 (0.1%)	Tablets
GOOD :	1656431 (99.9%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	pb97604	Batch NO. :	KQM2510440
Print Interval :	50000 Tablets	Running Time :	08 Hour 56 Min 18 Sec
Total Counter :	1700856	Total Running Time :	2999 Hour 24 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	92	%	Manual ;	45	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	36.99	kN	No		
HEP :	25.0 %	35.57	kN	Yes	10	Rotations 15 Tablets
HCP :	2.0 %	29.03	kN	Yes	15	Times
Ref :		28.46	kN			
LCP :	2.0 %	27.89	kN	Yes	15	Times
LEP :	25.0 %	21.34	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	19.92	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	3.95	kN	Min :	3.29	kN	Punch NO.;	23
SD :	8.50	%	Max :	4.75	kN	Punch NO.;	40
Mean Main Pressure :	27.52	kN	Min :	23.79	kN	Punch NO.;	23
SD :	6.90	%	Max :	31.79	kN	Punch NO.;	42
Filling Depth :	Increase ;	325	Times	Decrease ;	296	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	24	RPM	Capacity :	221760	Tablets/hour
Feeder : Status ;	Auto		Speed ;	22	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.20	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.85	mm	Main Thickness ;	0.85	mm
Filling Depth :	3.81	mm	Current Cam :	Cam B	
Main Air pressure :	0.642	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	45	Min			
Lower Punch Remain Time (S2):	7	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.62	kN		
	Max ;	0.85	kN	Max Punch ;	8
Lower Punch :	Average ;	0.45	kN		
	Max ;	1.10	kN	Max Punch ;	13
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	1700856	Tablets			
A.W.C Counter :	1657793	Tablets			
HEP :	334 (0.0%)	Tablets	LEP :	1028 (0.1%)	Tablets
GOOD :	1656431 (99.9%)	Tablets			