

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	sj97904	Batch NO. :	KQM2510440
Print Interval :	50000 Tablets	Running Time :	12 Hour 33 Min 55 Sec
Total Counter :	2501330	Total Running Time :	3003 Hour 2 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	82	%	Manual ;	45	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	41.55	kN	No		
HEP :	25.0 %	39.95	kN	Yes	10 Rotations	15 Tablets
HCP :	2.0 %	32.60	kN	Yes	15 Times	
Ref :		31.96	kN			
LCP :	2.0 %	31.32	kN	Yes	15 Times	
LEP :	25.0 %	23.97	kN	Yes	10 Rotations	15 Tablets
LSP :	30.0 %	22.37	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	4.17	kN	Min :	3.69	kN	Punch NO.;	42
SD :	6.20	%	Max :	4.79	kN	Punch NO.;	6
Mean Main Pressure :	32.36	kN	Min :	29.22	kN	Punch NO.;	42
SD :	4.90	%	Max :	35.38	kN	Punch NO.;	50
Filling Depth :	Increase ;	744	Times	Decrease ;	689	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	24	RPM	Capacity :	221760	Tablets/hour
Feeder : Status ;	Auto		Speed ;	20	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.20	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.85	mm	Main Thickness ;	0.85	mm
Filling Depth :	3.65	mm	Current Cam :	Cam B	
Main Air pressure :	0.750	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	28	Min			
Lower Punch Remain Time (S2):	1	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.41	kN		
	Max ;	0.49	kN	Max Punch ;	19
Lower Punch :	Average ;	0.42	kN		
	Max ;	1.27	kN	Max Punch ;	40
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	2501330	Tablets			
A.W.C Counter :	2452257	Tablets			
HEP :	511 (0.0%)	Tablets	LEP :	1393 (0.1%)	Tablets
GOOD :	2450353 (99.9%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	sj97904	Batch NO. :	KQM2510440
Print Interval :	50000 Tablets	Running Time :	12 Hour 33 Min 55 Sec
Total Counter :	2501330	Total Running Time :	3003 Hour 2 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	54	%	Manual ;	45	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	37.18	kN	No		
HEP :	25.0 %	35.75	kN	Yes	10	Rotations 15 Tablets
HCP :	2.0 %	29.17	kN	Yes	15	Times
Ref :		28.60	kN			
LCP :	2.0 %	28.03	kN	Yes	15	Times
LEP :	25.0 %	21.45	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	20.02	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	4.06	kN	Min :	3.56	kN	Punch NO.;	17
SD :	5.20	%	Max :	4.52	kN	Punch NO.;	32
Mean Main Pressure :	28.69	kN	Min :	25.77	kN	Punch NO.;	18
SD :	4.30	%	Max :	31.33	kN	Punch NO.;	32
Filling Depth :	Increase ;	561	Times	Decrease ;	533	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	24	RPM	Capacity :	221760	Tablets/hour
Feeder : Status ;	Auto		Speed ;	13	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.20	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.85	mm	Main Thickness ;	0.85	mm
Filling Depth :	3.70	mm	Current Cam :	Cam B 4.5-10mm	
Main Air pressure :	0.750	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	28	Min			
Lower Punch Remain Time (S2):	1	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.41	kN		
	Max ;	0.52	kN	Max Punch ;	26
Lower Punch :	Average ;	0.42	kN		
	Max ;	1.27	kN	Max Punch ;	40
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	2501330	Tablets			
A.W.C Counter :	2452257	Tablets			
HEP :	511 (0.0%)	Tablets	LEP :	1393 (0.1%)	Tablets
GOOD :	2450353 (99.9%)	Tablets			