

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	pb97604	Batch NO. :	KQM2510540
Print Interval :	50000 Tablets	Running Time :	00 Hour 52 Min 13 Sec
Total Counter :	200212	Total Running Time :	3004 Hour 38 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	93	%	Manual ;	45	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	40.77	kN	No		
HEP :	25.0 %	39.20	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	31.99	kN	Yes 15	Times	
Ref :		31.36	kN			
LCP :	2.0 %	30.73	kN	Yes 15	Times	
LEP :	25.0 %	23.52	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	21.95	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	3.92	kN	Min :	3.46	kN	Punch NO.;	39
SD :	5.40	%	Max :	4.43	kN	Punch NO.;	66
Mean Main Pressure :	31.50	kN	Min :	27.90	kN	Punch NO.;	39
SD :	4.80	%	Max :	35.35	kN	Punch NO.;	5
Filling Depth :	Increase ;	892	Times	Decrease ;	836	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	231000	Tablets/hour
Feeder : Status ;	Auto		Speed ;	23	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.20	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.85	mm	Main Thickness ;	0.85	mm
Filling Depth :	4.10	mm	Current Cam :	Cam B	
Main Air pressure :	0.710	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	35	Min			
Lower Punch Remain Time (S2):	10	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.58	kN		
	Max ;	0.72	kN	Max Punch ;	35
Lower Punch :	Average ;	0.29	kN		
	Max ;	0.88	kN	Max Punch ;	40
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	200212	Tablets			
A.W.C Counter :	192152	Tablets			
HEP :	1 (0.0%)	Tablets	LEP :	2 (0.0%)	Tablets
GOOD :	192149	Tablets			
	(100.0%)				

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	pb97604	Batch NO. :	KQM2510540
Print Interval :	50000 Tablets	Running Time :	00 Hour 52 Min 13 Sec
Total Counter :	200212	Total Running Time :	3004 Hour 38 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	87	%	Manual ;	45	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	37.08	kN	No		
HEP :	25.0 %	35.66	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	29.10	kN	Yes 15	Times	
Ref :		28.52	kN			
LCP :	2.0 %	27.95	kN	Yes 15	Times	
LEP :	25.0 %	21.39	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	19.97	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	3.86	kN	Min :	3.42	kN	Punch NO.;	77
SD :	4.60	%	Max :	4.27	kN	Punch NO.;	19
Mean Main Pressure :	28.46	kN	Min :	25.03	kN	Punch NO.;	77
SD :	4.60	%	Max :	31.13	kN	Punch NO.;	19
Filling Depth :	Increase ;	673	Times	Decrease ;	642	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	231000	Tablets/hour
Feeder : Status ;	Auto		Speed ;	22	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.20	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.85	mm	Main Thickness ;	0.85	mm
Filling Depth :	3.81	mm	Current Cam :	Cam B	
Main Air pressure :	0.710	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	35	Min			
Lower Punch Remain Time (S2):	10	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.58	kN		
	Max ;	0.66	kN	Max Punch ;	19
Lower Punch :	Average ;	0.29	kN		
	Max ;	0.88	kN	Max Punch ;	40
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	200212	Tablets			
A.W.C Counter :	192152	Tablets			
HEP :	1 (0.0%)	Tablets	LEP :	2 (0.0%)	Tablets
GOOD :	192149	Tablets			
	(100.0%)				