

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | METHO 500MG<br>STSH1A00 |
| User ID :        | ms88737       | Batch NO. :          | HAM2501350              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 34 Min 50 Sec   |
| Total Counter :  | 150064        | Total Running Time : | 3054 Hour 55 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 100       | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |        | 8.5-14mm  |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |        | 750000    | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 50        | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20        | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |        | 60.0      | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 3.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.90      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 23.0 % | 12.28     | kN      | No                 |              |           |
| HEP :                              | 18.0 % | 11.78     | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 1.5 %  | 10.13     | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 9.98      | kN      |                    |              |           |
| LCP :                              | 1.5 %  | 9.83      | kN      | Yes 15             | Times        |           |
| LEP :                              | 18.0 % | 8.19      | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 23.0 % | 7.69      | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 40.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None   |           |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.11       | kN  | Min : | 2.09       | kN  | Punch NO.; | 3  |
| SD :                 | 0.30       | %   | Max : | 2.12       | kN  | Punch NO.; | 34 |
| Mean Main Pressure : | 9.92       | kN  | Min : | 9.08       | kN  | Punch NO.; | 2  |
| SD :                 | 3.50       | %   | Max : | 10.63      | kN  | Punch NO.; | 9  |
| Filling Depth :      | Increase ; | 128 | Times | Decrease ; | 141 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 15       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 8.00    | mm  | Pre Thickness ;  | 5.00     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.65    | mm  | Main Thickness ; | 2.65     | mm           |
| Filling Depth :               | 8.48    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.740   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 2       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 26 |
| Lower Punch :    | Average ; | 0.13 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 12 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 150064            | Tablets |       |            |         |
| A.W.C Counter : | 129257            | Tablets |       |            |         |
| HEP :           | 815 (0.6%)        | Tablets | LEP : | 882 (0.7%) | Tablets |
| GOOD :          | 127560<br>(98.7%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | METHO 500MG<br>STSH1A00 |
| User ID :        | ms88737       | Batch NO. :          | HAM2501350              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 34 Min 50 Sec   |
| Total Counter :  | 150064        | Total Running Time : | 3054 Hour 55 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |           |         |                    |       |                      |
|------------------------------------|--------|-----------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 100       | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min : | 6.50 mm              |
|                                    |        | 8.5-14mm  |         |                    | Max : | 12.00 mm             |
| Target Quantity :                  |        | 750000    | Tablets |                    |       |                      |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |           |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 50        | Min     | Supply Time ;      | 0.2   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 20        | Min     | Supply Time ;      | 3.0   | Sec                  |
| Powder Supply Time :               |        | 60.0      | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 1         | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 2.00      | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 3.00      | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.90      | kN      |                    |       |                      |
|                                    |        |           |         | Stop Condition ;   |       |                      |
| HSP :                              | 23.0 % | 11.20     | kN      | No                 |       |                      |
| HEP :                              | 18.0 % | 10.74     | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 1.5 %  | 9.24      | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 9.10      | kN      |                    |       |                      |
| LCP :                              | 1.5 %  | 8.97      | kN      | Yes                | 15    | Times                |
| LEP :                              | 18.0 % | 7.47      | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 23.0 % | 7.01      | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |           |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00     | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |           |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.13       | kN  | Min : | 2.11       | kN  | Punch NO.; | 33 |
| SD :                 | 0.40       | %   | Max : | 2.15       | kN  | Punch NO.; | 18 |
| Mean Main Pressure : | 9.13       | kN  | Min : | 8.49       | kN  | Punch NO.; | 28 |
| SD :                 | 3.70       | %   | Max : | 9.78       | kN  | Punch NO.; | 22 |
| Filling Depth :      | Increase ; | 143 | Times | Decrease ; | 167 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 15       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 8.00    | mm  | Pre Thickness ;  | 5.00     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.85    | mm  | Main Thickness ; | 2.85     | mm           |
| Filling Depth :               | 8.68    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.740   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 2       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 26 |
| Lower Punch :    | Average ; | 0.13 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 12 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 150064            | Tablets |       |            |         |
| A.W.C Counter : | 129257            | Tablets |       |            |         |
| HEP :           | 815 (0.6%)        | Tablets | LEP : | 882 (0.7%) | Tablets |
| GOOD :          | 127560<br>(98.7%) | Tablets |       |            |         |