

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICL 200MG,STVI1A00  |
| User ID :        | rp98927       | Batch NO. :          | VIM2500620            |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 21 Min 03 Sec |
| Total Counter :  | 50028         | Total Running Time : | 3084 Hour 45 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 90       | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 1740000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 50       | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 23.0 %      | 11.20    | kN      | No                 |              |           |
| HEP :                              | 20.0 %      | 10.93    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 1.5 %       | 9.24     | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 9.10     | kN      |                    |              |           |
| LCP :                              | 1.5 %       | 8.97     | kN      | Yes 15             | Times        |           |
| LEP :                              | 20.0 %      | 7.28     | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 23.0 %      | 7.01     | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |              |       |            |       |       |            |    |
|----------------------|--------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 2.32         | kN    | Min :      | 2.30  | kN    | Punch NO.; | 14 |
| SD :                 | 0.40         | %     | Max :      | 2.35  | kN    | Punch NO.; | 6  |
| Mean Main Pressure : | 9.04         | kN    | Min :      | 8.33  | kN    | Punch NO.; | 9  |
| SD :                 | 4.00         | %     | Max :      | 10.17 | kN    | Punch NO.; | 4  |
| Filling Depth :      | Increase ; 3 | Times | Decrease ; | 13    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 156000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 18       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.50    | mm  | Pre Thickness ;  | 3.51     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.26    | mm  | Main Thickness ; | 2.27     | mm           |
| Filling Depth :               | 6.10    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.781   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 40      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 13      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.19 | kN |             |    |
|                  | Max ;     | 0.22 | kN | Max Punch ; | 52 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 40 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                  |         |       |        |         |
|-----------------|------------------|---------|-------|--------|---------|
| Total Counter : | 50028            | Tablets |       |        |         |
| A.W.C Counter : | 43361            | Tablets |       |        |         |
| HEP :           | 33 (0.1%)        | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 43328<br>(99.9%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICL 200MG,STVI1A00  |
| User ID :        | rp98927       | Batch NO. :          | VIM2500620            |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 21 Min 03 Sec |
| Total Counter :  | 50028         | Total Running Time : | 3084 Hour 45 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 90        | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |        | Cam B     |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |        | 4.5-10mm  |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |        | 1740000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 50        | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20        | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |        | 60.0      | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 3.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.90      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 23.0 % | 11.98     | kN      | No                 |              |           |
| HEP :                              | 20.0 % | 11.69     | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 1.5 %  | 9.88      | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 9.74      | kN      |                    |              |           |
| LCP :                              | 1.5 %  | 9.59      | kN      | Yes 15             | Times        |           |
| LEP :                              | 20.0 % | 7.79      | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 23.0 % | 7.50      | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 40.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |           |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.35       | kN | Min : | 2.32       | kN | Punch NO.; | 24 |
| SD :                 | 0.50       | %  | Max : | 2.38       | kN | Punch NO.; | 39 |
| Mean Main Pressure : | 9.76       | kN | Min : | 9.31       | kN | Punch NO.; | 12 |
| SD :                 | 3.00       | %  | Max : | 10.77      | kN | Punch NO.; | 20 |
| Filling Depth :      | Increase ; | 13 | Times | Decrease ; | 12 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 156000            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 18                | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 6.50    | mm  | Pre Thickness ;  | 3.51              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.45    | mm  | Main Thickness ; | 2.46              | mm           |
| Filling Depth :               | 6.00    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.781   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 40      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 13      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.19 | kN |             |    |
|                  | Max ;     | 0.21 | kN | Max Punch ; | 52 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 40 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                  |         |       |        |         |
|-----------------|------------------|---------|-------|--------|---------|
| Total Counter : | 50028            | Tablets |       |        |         |
| A.W.C Counter : | 43361            | Tablets |       |        |         |
| HEP :           | 33 (0.1%)        | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 43328<br>(99.9%) | Tablets |       |        |         |