

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICL 200MG,STVI1A00  |
| User ID :        | rp98927       | Batch NO. :          | VIM2500620            |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 16 Min 48 Sec |
| Total Counter :  | 350170        | Total Running Time : | 3086 Hour 40 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 90       | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 1740000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 50       | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 23.0 %      | 11.20    | kN      | No                 |              |           |
| HEP :                              | 20.0 %      | 10.93    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 1.5 %       | 9.24     | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 9.10     | kN      |                    |              |           |
| LCP :                              | 1.5 %       | 8.97     | kN      | Yes 15             | Times        |           |
| LEP :                              | 20.0 %      | 7.28     | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 23.0 %      | 7.01     | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |              |       |            |       |       |            |    |
|----------------------|--------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 2.38         | kN    | Min :      | 2.34  | kN    | Punch NO.; | 1  |
| SD :                 | 0.70         | %     | Max :      | 2.43  | kN    | Punch NO.; | 3  |
| Mean Main Pressure : | 9.12         | kN    | Min :      | 8.53  | kN    | Punch NO.; | 64 |
| SD :                 | 4.20         | %     | Max :      | 10.47 | kN    | Punch NO.; | 3  |
| Filling Depth :      | Increase ; 9 | Times | Decrease ; | 21    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 156000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 18       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.30    | mm  | Pre Thickness ;  | 3.31     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.26    | mm  | Main Thickness ; | 2.27     | mm           |
| Filling Depth :               | 6.04    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.780   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 25      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 18      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.20 | kN |             |    |
|                  | Max ;     | 0.23 | kN | Max Punch ; | 1  |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 21 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 350170             | Tablets |       |        |         |
| A.W.C Counter : | 343055             | Tablets |       |        |         |
| HEP :           | 39 (0.0%)          | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 343016<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICL 200MG,STVI1A00  |
| User ID :        | rp98927       | Batch NO. :          | VIM2500620            |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 16 Min 48 Sec |
| Total Counter :  | 350170        | Total Running Time : | 3086 Hour 40 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 90        | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |        | Cam B     |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |        | 4.5-10mm  |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |        | 1740000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 50        | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20        | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |        | 60.0      | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 3.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.90      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 23.0 % | 11.40     | kN      | No                 |              |           |
| HEP :                              | 20.0 % | 11.12     | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 1.5 %  | 9.41      | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 9.27      | kN      |                    |              |           |
| LCP :                              | 1.5 %  | 9.13      | kN      | Yes 15             | Times        |           |
| LEP :                              | 20.0 % | 7.42      | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 23.0 % | 7.14      | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 40.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |           |         |                    |              |           |

### Pressure Data

|                      |               |       |               |       |    |            |    |
|----------------------|---------------|-------|---------------|-------|----|------------|----|
| Mean Pre Pressure :  | 2.48          | kN    | Min :         | 2.45  | kN | Punch NO.; | 57 |
| SD :                 | 0.60          | %     | Max :         | 2.52  | kN | Punch NO.; | 4  |
| Mean Main Pressure : | 9.29          | kN    | Min :         | 8.79  | kN | Punch NO.; | 52 |
| SD :                 | 3.10          | %     | Max :         | 10.07 | kN | Punch NO.; | 6  |
| Filling Depth :      | Increase ; 27 | Times | Decrease ; 26 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 156000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 18       | RPM          |
| Pre Upper Punch Position ;    | 2.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.20    | mm  | Pre Thickness ;  | 3.20     | mm           |
| Main Upper Punch Position ;   | 2.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.50    | mm  | Main Thickness ; | 2.50     | mm           |
| Filling Depth :               | 6.00    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.780   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 25      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 18      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.20 | kN |             |    |
|                  | Max ;     | 0.23 | kN | Max Punch ; | 64 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 21 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 350170             | Tablets |       |        |         |
| A.W.C Counter : | 343055             | Tablets |       |        |         |
| HEP :           | 39 (0.0%)          | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 343016<br>(100.0%) | Tablets |       |        |         |