

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Allo300mg STUB2B00    |
| User ID :        | rp98927       | Batch NO. :          | UBM2501830            |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 45 Min 54 Sec |
| Total Counter :  | 350080        | Total Running Time : | 3099 Hour 6 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 71       | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 1000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 50       | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 23.0 %      | 12.42    | kN      | No                 |              |           |
| HEP :                              | 20.0 %      | 12.12    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 1.5 %       | 10.25    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 10.10    | kN      |                    |              |           |
| LCP :                              | 1.5 %       | 9.95     | kN      | Yes 15             | Times        |           |
| LEP :                              | 20.0 %      | 8.08     | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 23.0 %      | 7.78     | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 1.01           | kN    | Min :      | 0.90  | kN    | Punch NO.; | 29 |
| SD :                 | 4.40           | %     | Max :      | 1.11  | kN    | Punch NO.; | 19 |
| Mean Main Pressure : | 10.15          | kN    | Min :      | 9.04  | kN    | Punch NO.; | 42 |
| SD :                 | 4.80           | %     | Max :      | 11.26 | kN    | Punch NO.; | 52 |
| Filling Depth :      | Increase ; 254 | Times | Decrease ; | 264   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 14       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.00    | mm  | Pre Thickness ;  | 4.00     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.15    | mm  | Main Thickness ; | 3.15     | mm           |
| Filling Depth :               | 6.80    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.780   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 35      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 19      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 18 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |          |         |
|-----------------|--------------------|---------|-------|----------|---------|
| Total Counter : | 350080             | Tablets |       |          |         |
| A.W.C Counter : | 345925             | Tablets |       |          |         |
| HEP :           | 25 (0.0%)          | Tablets | LEP : | 1 (0.0%) | Tablets |
| GOOD :          | 345899<br>(100.0%) | Tablets |       |          |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Allo300mg STUB2B00    |
| User ID :        | rp98927       | Batch NO. :          | UBM2501830            |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 45 Min 54 Sec |
| Total Counter :  | 350080        | Total Running Time : | 3099 Hour 6 Min       |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 65       | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |             | 1000000  | Tablets |                    |              |           |
| Oil Lubrication :                  |             | 0.000    | MPa     |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 50       | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 23.0 %      | 11.40    | kN      | No                 |              |           |
| HEP :                              | 20.0 %      | 11.12    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 1.5 %       | 9.41     | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 9.27     | kN      |                    |              |           |
| LCP :                              | 1.5 %       | 9.13     | kN      | Yes 15             | Times        |           |
| LEP :                              | 20.0 %      | 7.42     | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 23.0 %      | 7.14     | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 1.08           | kN    | Min :      | 0.94  | kN    | Punch NO.; | 32 |
| SD :                 | 4.80           | %     | Max :      | 1.19  | kN    | Punch NO.; | 15 |
| Mean Main Pressure : | 9.37           | kN    | Min :      | 7.74  | kN    | Punch NO.; | 42 |
| SD :                 | 5.50           | %     | Max :      | 10.41 | kN    | Punch NO.; | 52 |
| Filling Depth :      | Increase ; 234 | Times | Decrease ; | 218   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 13       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.00    | mm  | Pre Thickness ;  | 4.00     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.30    | mm  | Main Thickness ; | 3.30     | mm           |
| Filling Depth :               | 6.90    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.780   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 35      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 19      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 18 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |          |         |
|-----------------|--------------------|---------|-------|----------|---------|
| Total Counter : | 350080             | Tablets |       |          |         |
| A.W.C Counter : | 345925             | Tablets |       |          |         |
| HEP :           | 25 (0.0%)          | Tablets | LEP : | 1 (0.0%) | Tablets |
| GOOD :          | 345899<br>(100.0%) | Tablets |       |          |         |