

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | METHO<br>750MG,STHA1B00 |
| User ID :        | rp98927       | Batch NO. :          | HAM2500975              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 34 Min 58 Sec   |
| Total Counter :  | 150050        | Total Running Time : | 3121 Hour 25 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 201               | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min : | 6.50 mm              |
| Target Quantity :                  |        | 500000            | Tablets |                    | Max : | 12.00 mm             |
| Air Pressure Low Limit :           |        | 0.000             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 30                | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 3.0   | Sec                  |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 1                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 3.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.90              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 11.00             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 10.57             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 1.5 %  | 8.58              | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 8.46              | kN      |                    |       |                      |
| LCP :                              | 1.5 %  | 8.33              | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 6.34              | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 5.92              | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |      |       |            |      |            |    |
|----------------------|------------|------|-------|------------|------|------------|----|
| Mean Pre Pressure :  | 1.56       | kN   | Min : | 1.44       | kN   | Punch NO.; | 42 |
| SD :                 | 3.20       | %    | Max : | 1.73       | kN   | Punch NO.; | 17 |
| Mean Main Pressure : | 8.66       | kN   | Min : | 7.94       | kN   | Punch NO.; | 40 |
| SD :                 | 2.90       | %    | Max : | 9.48       | kN   | Punch NO.; | 17 |
| Filling Depth :      | Increase ; | 1537 | Times | Decrease ; | 1526 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400             | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 30                | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 8.20    | mm  | Pre Thickness ;  | 5.20              | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 6.95    | mm  | Main Thickness ; | 3.95              | mm           |
| Filling Depth :               | 9.56    | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.748   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 8       | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.23 | kN |             |    |
|                  | Max ;     | 0.28 | kN | Max Punch ; | 29 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 4  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 150050            | Tablets |       |            |         |
| A.W.C Counter : | 141761            | Tablets |       |            |         |
| HEP :           | 191 (0.1%)        | Tablets | LEP : | 124 (0.1%) | Tablets |
| GOOD :          | 141446<br>(99.8%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | METHO<br>750MG,STHA1B00 |
| User ID :        | rp98927       | Batch NO. :          | HAM2500975              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 34 Min 58 Sec   |
| Total Counter :  | 150050        | Total Running Time : | 3121 Hour 25 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 193               | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |        | 500000            | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |        | 0.000             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 30                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.90              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 11.61             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 11.17             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 1.5 %  | 9.07              | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 8.93              | kN      |                    |              |           |
| LCP :                              | 1.5 %  | 8.80              | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 6.70              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 6.25              | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |      |       |            |      |            |    |
|----------------------|------------|------|-------|------------|------|------------|----|
| Mean Pre Pressure :  | 1.28       | kN   | Min : | 1.18       | kN   | Punch NO.; | 42 |
| SD :                 | 4.70       | %    | Max : | 1.42       | kN   | Punch NO.; | 19 |
| Mean Main Pressure : | 9.04       | kN   | Min : | 8.52       | kN   | Punch NO.; | 10 |
| SD :                 | 3.70       | %    | Max : | 10.14      | kN   | Punch NO.; | 19 |
| Filling Depth :      | Increase ; | 1215 | Times | Decrease ; | 1184 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400             | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 29                | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 8.05    | mm  | Pre Thickness ;  | 5.05              | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 7.15    | mm  | Main Thickness ; | 4.15              | mm           |
| Filling Depth :               | 9.65    | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.748   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 8       | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.23 | kN |             |    |
|                  | Max ;     | 0.28 | kN | Max Punch ; | 29 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 4  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 150050            | Tablets |       |            |         |
| A.W.C Counter : | 141761            | Tablets |       |            |         |
| HEP :           | 191 (0.1%)        | Tablets | LEP : | 124 (0.1%) | Tablets |
| GOOD :          | 141446<br>(99.8%) | Tablets |       |            |         |