

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA50mg<br>STSH4A00 |
| User ID :        | pb97604       | Batch NO. :          | SHM2503850             |
| Print Interval : | 50000 Tablets | Running Time :       | 12 Hour 11 Min 48 Sec  |
| Total Counter :  | 3102340       | Total Running Time : | 3152 Hour 14 Min       |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 28                | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 6000000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 30                | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 30                | Min     | Supply Time ;      | 3.0   | Sec                  |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 1                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 3.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.90              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 17.67             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 16.99             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 1.5 %  | 13.80             | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 13.59             | kN      |                    |       |                      |
| LCP :                              | 1.5 %  | 13.39             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 10.20             | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 9.52              | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |     |            |    |
|----------------------|------------|----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.54       | kN | Min : | 1.51       | kN  | Punch NO.; | 33 |
| SD :                 | 0.80       | %  | Max : | 1.57       | kN  | Punch NO.; | 5  |
| Mean Main Pressure : | 13.72      | kN | Min : | 12.44      | kN  | Punch NO.; | 25 |
| SD :                 | 3.90       | %  | Max : | 15.42      | kN  | Punch NO.; | 39 |
| Filling Depth :      | Increase ; | 65 | Times | Decrease ; | 131 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 8        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.70    | mm  | Pre Thickness ;  | 2.71     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.06    | mm  | Main Thickness ; | 2.07     | mm           |
| Filling Depth :               | 5.44    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.746   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 16      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.18 | kN |             |    |
|                  | Max ;     | 0.24 | kN | Max Punch ; | 70 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 69 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 3102340             | Tablets |       |        |         |
| A.W.C Counter : | 3096568             | Tablets |       |        |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 3096568<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA50mg<br>STSH4A00 |
| User ID :        | pb97604       | Batch NO. :          | SHM2503850             |
| Print Interval : | 50000 Tablets | Running Time :       | 12 Hour 11 Min 48 Sec  |
| Total Counter :  | 3102340       | Total Running Time : | 3152 Hour 14 Min       |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 27                | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 6000000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 30                | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 30                | Min     | Supply Time ;      | 3.0   | Sec                  |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 1                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 3.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.90              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 18.24             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 17.53             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 1.5 %  | 14.24             | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 14.03             | kN      |                    |       |                      |
| LCP :                              | 1.5 %  | 13.82             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 10.52             | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 9.82              | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1     | Rotations            |
| Empty Punch No. :                  | N      |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 1.46       | kN | Min : | 1.43       | kN | Punch NO.; | 59 |
| SD :                 | 0.90       | %  | Max : | 1.50       | kN | Punch NO.; | 21 |
| Mean Main Pressure : | 14.00      | kN | Min : | 12.82      | kN | Punch NO.; | 74 |
| SD :                 | 4.20       | %  | Max : | 15.85      | kN | Punch NO.; | 60 |
| Filling Depth :      | Increase ; | 34 | Times | Decrease ; | 68 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 8        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.70    | mm  | Pre Thickness ;  | 2.71     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.15    | mm  | Main Thickness ; | 2.16     | mm           |
| Filling Depth :               | 5.54    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.746   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 16      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.18 | kN |             |    |
|                  | Max ;     | 0.25 | kN | Max Punch ; | 62 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 69 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 3102340             | Tablets |       |        |         |
| A.W.C Counter : | 3096568             | Tablets |       |        |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 3096568<br>(100.0%) | Tablets |       |        |         |