

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA50mg<br>STSH4A00 |
| User ID :        | pb97604       | Batch NO. :          | SHM2503950             |
| Print Interval : | 50000 Tablets | Running Time :       | 05 Hour 24 Min 00 Sec  |
| Total Counter :  | 1401212       | Total Running Time : | 3168 Hour 17 Min       |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |       |                      |
|------------------------------------|-------------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ;      | 28                | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |             | 6000000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |             |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |             | 30                | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |             | 30                | Min     | Supply Time ;      | 3.0   | Sec                  |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |       |                      |
| Initial Reject Time :              |             | 1                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |             | 3.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |             | 4.90              | kN      |                    |       |                      |
|                                    |             |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 %      | 18.19             | kN      | No                 |       |                      |
| HEP :                              | 25.0 %      | 17.49             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 1.5 %       | 14.20             | kN      | Yes                | 15    | Times                |
| Ref :                              |             | 13.99             | kN      |                    |       |                      |
| LCP :                              | 1.5 %       | 13.78             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 %      | 10.50             | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 %      | 9.80              | kN      | No                 |       |                      |
| SD limit :                         | 45.0 %      |                   |         | Yes                |       |                      |
| Pre HSP                            |             | 40.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 1     | Rotations            |
| Empty Punch No. :                  | None        |                   |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.56       | kN  | Min : | 1.52       | kN  | Punch NO.; | 42 |
| SD :                 | 0.70       | %   | Max : | 1.59       | kN  | Punch NO.; | 37 |
| Mean Main Pressure : | 13.95      | kN  | Min : | 12.58      | kN  | Punch NO.; | 74 |
| SD :                 | 3.80       | %   | Max : | 15.13      | kN  | Punch NO.; | 53 |
| Filling Depth :      | Increase ; | 133 | Times | Decrease ; | 199 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 8        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.70    | mm  | Pre Thickness ;  | 2.71     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.06    | mm  | Main Thickness ; | 2.07     | mm           |
| Filling Depth :               | 5.37    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.768   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 18      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 8       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.20 | kN |             |    |
|                  | Max ;     | 0.27 | kN | Max Punch ; | 4  |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 25 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 1401212             | Tablets |       |        |         |
| A.W.C Counter : | 1396413             | Tablets |       |        |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 1396413<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA50mg<br>STSH4A00 |
| User ID :        | pb97604       | Batch NO. :          | SHM2503950             |
| Print Interval : | 50000 Tablets | Running Time :       | 05 Hour 24 Min 00 Sec  |
| Total Counter :  | 1401212       | Total Running Time : | 3168 Hour 17 Min       |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 27                | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 6000000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 30                | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 30                | Min     | Supply Time ;      | 3.0   | Sec                  |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 1                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 3.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.90              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 18.40             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 17.70             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 1.5 %  | 14.37             | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 14.16             | kN      |                    |       |                      |
| LCP :                              | 1.5 %  | 13.94             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 10.62             | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 9.91              | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1     | Rotations            |
| Empty Punch No. :                  | N      |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 1.47       | kN | Min : | 1.43       | kN | Punch NO.; | 74 |
| SD :                 | 1.00       | %  | Max : | 1.51       | kN | Punch NO.; | 55 |
| Mean Main Pressure : | 14.17      | kN | Min : | 12.30      | kN | Punch NO.; | 6  |
| SD :                 | 4.30       | %  | Max : | 15.84      | kN | Punch NO.; | 60 |
| Filling Depth :      | Increase ; | 59 | Times | Decrease ; | 93 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 8        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.70    | mm  | Pre Thickness ;  | 2.71     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.15    | mm  | Main Thickness ; | 2.16     | mm           |
| Filling Depth :               | 5.49    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.768   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 18      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 8       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.20 | kN |             |    |
|                  | Max ;     | 0.27 | kN | Max Punch ; | 4  |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 19 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 1401212             | Tablets |       |        |         |
| A.W.C Counter : | 1396413             | Tablets |       |        |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 1396413<br>(100.0%) | Tablets |       |        |         |