

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                          |
|------------------|---------------|----------------------|--------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRALINE<br>100mg TABS |
| User ID :        | vg96371       | Batch NO. :          | SED-TAB-XV-25-001B       |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 23 Min 37 Sec    |
| Total Counter :  | 550258        | Total Running Time : | 3190 Hour 58 Min         |
| Station :        | Station 1     |                      |                          |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 30                | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 650000            | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 30                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30                | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 14.89             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 14.32             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 1.5 %       | 11.63             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 11.46             | kN      |                    |              |           |
| LCP :                              | 1.5 %       | 11.28             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 8.59              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 8.02              | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 1.44           | kN    | Min :      | 1.42  | kN    | Punch NO.; | 68 |
| SD :                 | 0.70           | %     | Max :      | 1.47  | kN    | Punch NO.; | 2  |
| Mean Main Pressure : | 11.30          | kN    | Min :      | 10.78 | kN    | Punch NO.; | 18 |
| SD :                 | 1.70           | %     | Max :      | 11.65 | kN    | Punch NO.; | 24 |
| Filling Depth :      | Increase ; 355 | Times | Decrease ; | 316   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 231000            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 8                 | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 8.00    | mm  | Pre Thickness ;  | 5.01              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 6.08    | mm  | Main Thickness ; | 3.09              | mm           |
| Filling Depth :               | 7.99    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.762   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 21      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 12      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.23 | kN | Max Punch ; | 1  |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 74 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |           |         |
|-----------------|--------------------|---------|-------|-----------|---------|
| Total Counter : | 550258             | Tablets |       |           |         |
| A.W.C Counter : | 538924             | Tablets |       |           |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 42 (0.0%) | Tablets |
| GOOD :          | 538882<br>(100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                          |
|------------------|---------------|----------------------|--------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRALINE<br>100mg TABS |
| User ID :        | vg96371       | Batch NO. :          | SED-TAB-XV-25-001B       |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 23 Min 37 Sec    |
| Total Counter :  | 550258        | Total Running Time : | 3190 Hour 58 Min         |
| Station :        | Station 2     |                      |                          |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 30                | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 650000            | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 30                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 30                | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.90              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 15.89             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 15.28             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 1.5 %  | 12.41             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 12.23             | kN      |                    |              |           |
| LCP :                              | 1.5 %  | 12.04             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 9.17              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 8.56              | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.35  | kN    | Min :      | 1.29  | kN    | Punch NO.; | 74 |
|                      |            | 1.80  | %     | Max :      | 1.40  | kN    | Punch NO.; | 17 |
| Mean Main Pressure : | SD :       | 12.29 | kN    | Min :      | 11.64 | kN    | Punch NO.; | 2  |
|                      |            | 2.80  | %     | Max :      | 13.61 | kN    | Punch NO.; | 66 |
| Filling Depth :      | Increase ; | 140   | Times | Decrease ; | 153   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 231000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 8        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 8.00    | mm  | Pre Thickness ;  | 5.01     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.23    | mm  | Main Thickness ; | 3.24     | mm           |
| Filling Depth :               | 7.68    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.762   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 21      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 12      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.23 | kN | Max Punch ; | 1  |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 74 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |           |         |
|-----------------|--------------------|---------|-------|-----------|---------|
| Total Counter : | 550258             | Tablets |       |           |         |
| A.W.C Counter : | 539001             | Tablets |       |           |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 42 (0.0%) | Tablets |
| GOOD :          | 538959<br>(100.0%) | Tablets |       |           |         |