

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ATORV 40<br>MG,STJB1C00 |
| User ID :        | gp98929       | Batch NO. :          | JBM2503640              |
| Print Interval : | 50000 Tablets | Running Time :       | 06 Hour 04 Min 18 Sec   |
| Total Counter :  | 750322        | Total Running Time : | 3249 Hour 47 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |                        |         |                    |       |                      |
|------------------------------------|--------|------------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 223                    | %       | Manual ;           | 27    | RPM                  |
| Filling Cam :                      |        | Cam D<br>12.5-<br>18mm |         | Cam Range :        | Min : | 10.50 mm             |
| Target Quantity :                  |        | 1500000                | Tablets |                    | Max : | 16.00 mm             |
| Air Pressure Low Limit :           |        | 0.400                  | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                        |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 50                     | Min     | Supply Time ;      | 0.2   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 30                     | Min     | Supply Time ;      | 2.0   | Sec                  |
| Powder Supply Time :               |        | 60.0                   | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                      | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 2.00                   | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 3.00                   | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.90                   | kN      |                    |       |                      |
|                                    |        |                        |         | Stop Condition ;   |       |                      |
| HSP :                              | 25.0 % | 12.12                  | kN      | No                 |       |                      |
| HEP :                              | 20.0 % | 11.63                  | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 9.89                   | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 9.69                   | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 9.50                   | kN      | Yes                | 15    | Times                |
| LEP :                              | 20.0 % | 7.76                   | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 25.0 % | 7.27                   | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |                        |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00                  | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations              |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |                        |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.84       | kN  | Min : | 1.76       | kN  | Punch NO.; | 46 |
| SD :                 | 1.50       | %   | Max : | 1.92       | kN  | Punch NO.; | 1  |
| Mean Main Pressure : | 10.61      | kN  | Min : | 9.01       | kN  | Punch NO.; | 27 |
| SD :                 | 10.40      | %   | Max : | 13.36      | kN  | Punch NO.; | 44 |
| Filling Depth :      | Increase ; | 318 | Times | Decrease ; | 383 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                        |              |
|-------------------------------|---------|-----|------------------|------------------------|--------------|
| Disk Speed :                  | 18      | RPM | Capacity :       | 114480                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 40                     | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                        |              |
| Pre Lower Punch Position ;    | 7.50    | mm  | Pre Thickness ;  | 4.50                   | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                        |              |
| Main Lower Punch Position ;   | 6.86    | mm  | Main Thickness ; | 3.86                   | mm           |
|                               |         |     |                  |                        |              |
| Filling Depth :               | 12.54   | mm  | Current Cam :    | Cam D<br>12.5-<br>18mm |              |
| Main Air pressure :           | 0.772   | MPa |                  |                        |              |
| Oil Lubrication status :      |         |     |                  |                        |              |
| Upper Punch Remain Time (S1): | 36      | Min |                  |                        |              |
| Lower Punch Remain Time (S2): | 13      | Min |                  |                        |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                        |              |
| Powder Status :               | Disable |     | Dust Collector : | Off                    |              |
| Initial Reject :              | On      |     | Buzzer :         | Off                    |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.01 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 49 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 47 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |           |         |
|-----------------|--------------------|---------|-------|-----------|---------|
| Total Counter : | 750322             | Tablets |       |           |         |
| A.W.C Counter : | 741451             | Tablets |       |           |         |
| HEP :           | 58 (0.0%)          | Tablets | LEP : | 44 (0.0%) | Tablets |
| GOOD :          | 741349<br>(100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ATORV 40<br>MG,STJB1C00 |
| User ID :        | gp98929       | Batch NO. :          | JBM2503640              |
| Print Interval : | 50000 Tablets | Running Time :       | 06 Hour 04 Min 18 Sec   |
| Total Counter :  | 750322        | Total Running Time : | 3249 Hour 47 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                        |         |                    |       |                      |
|------------------------------------|--------|------------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 223                    | %       | Manual ;           | 40    | RPM                  |
| Filling Cam :                      |        | Cam D<br>12.5-<br>18mm |         | Cam Range :        | Min : | 10.50 mm             |
| Target Quantity :                  |        | 1500000                | Tablets |                    | Max : | 16.00 mm             |
| Air Pressure Low Limit :           |        | 0.400                  | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                        |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 50                     | Min     | Supply Time ;      | 0.2   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 30                     | Min     | Supply Time ;      | 2.0   | Sec                  |
| Powder Supply Time :               |        | 60.0                   | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                      | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 2.00                   | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 3.00                   | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.90                   | kN      |                    |       |                      |
|                                    |        |                        |         | Stop Condition ;   |       |                      |
| HSP :                              | 25.0 % | 11.53                  | kN      | No                 |       |                      |
| HEP :                              | 20.0 % | 11.07                  | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 9.41                   | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 9.22                   | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 9.04                   | kN      | Yes                | 15    | Times                |
| LEP :                              | 20.0 % | 7.38                   | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 25.0 % | 6.92                   | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |                        |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00                  | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations              |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |                        |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.63       | kN  | Min : | 1.54       | kN  | Punch NO.; | 15 |
| SD :                 | 2.30       | %   | Max : | 1.73       | kN  | Punch NO.; | 22 |
| Mean Main Pressure : | 9.33       | kN  | Min : | 8.70       | kN  | Punch NO.; | 41 |
| SD :                 | 3.80       | %   | Max : | 10.09      | kN  | Punch NO.; | 5  |
| Filling Depth :      | Increase ; | 304 | Times | Decrease ; | 308 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                        |              |
|-------------------------------|---------|-----|------------------|------------------------|--------------|
| Disk Speed :                  | 18      | RPM | Capacity :       | 114480                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 40                     | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                        |              |
| Pre Lower Punch Position ;    | 7.50    | mm  | Pre Thickness ;  | 4.50                   | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                        |              |
| Main Lower Punch Position ;   | 7.05    | mm  | Main Thickness ; | 4.05                   | mm           |
|                               |         |     |                  |                        |              |
| Filling Depth :               | 12.31   | mm  | Current Cam :    | Cam D<br>12.5-<br>18mm |              |
| Main Air pressure :           | 0.772   | MPa |                  |                        |              |
| Oil Lubrication status :      |         |     |                  |                        |              |
| Upper Punch Remain Time (S1): | 36      | Min |                  |                        |              |
| Lower Punch Remain Time (S2): | 13      | Min |                  |                        |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                        |              |
| Powder Status :               | Disable |     | Dust Collector : | Off                    |              |
| Initial Reject :              | On      |     | Buzzer :         | Off                    |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.01 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 44 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 47 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |           |         |
|-----------------|--------------------|---------|-------|-----------|---------|
| Total Counter : | 750322             | Tablets |       |           |         |
| A.W.C Counter : | 741451             | Tablets |       |           |         |
| HEP :           | 58 (0.0%)          | Tablets | LEP : | 44 (0.0%) | Tablets |
| GOOD :          | 741349<br>(100.0%) | Tablets |       |           |         |