

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	GEMFI 600MG,STKF2A00
User ID :	rp98927	Batch NO. :	KFM2500260
Print Interval :	50000 Tablets	Running Time :	01 Hour 03 Min 29 Sec
Total Counter :	100008	Total Running Time :	3293 Hour 35 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	174	%	Manual ;	15	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		600000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	0.6	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	4.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	9.66	kN	No		
HEP :	25.0 %	9.29	kN	Yes	10	Rotations 15 Tablets
HCP :	2.0 %	7.58	kN	Yes	15	Times
Ref :		7.43	kN			
LCP :	2.0 %	7.28	kN	Yes	15	Times
LEP :	25.0 %	5.57	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	5.20	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	1.42	kN	Min :	1.29	kN	Punch NO.;	44
SD :	4.10	%	Max :	1.55	kN	Punch NO.;	5
Mean Main Pressure :	7.43	kN	Min :	6.97	kN	Punch NO.;	30
SD :	3.20	%	Max :	7.97	kN	Punch NO.;	6
Filling Depth :	Increase ;	85	Times	Decrease ;	78	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	26	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	6.90	mm	Pre Thickness ;	3.90	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	6.55	mm	Main Thickness ;	3.55	mm
Filling Depth :	9.09	mm	Current Cam :	Cam C	
Main Air pressure :	0.752	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	40	Min			
Lower Punch Remain Time (S2):	1	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0
Lower Punch :	Average ;	0.15	kN		
	Max ;	0.29	kN	Max Punch ;	20
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	100008	Tablets			
A.W.C Counter :	93744	Tablets			
HEP :	166 (0.2%)	Tablets	LEP :	36 (0.0%)	Tablets
GOOD :	93542 (99.8%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	GEMFI 600MG,STKF2A00
User ID :	rp98927	Batch NO. :	KFM2500260
Print Interval :	50000 Tablets	Running Time :	01 Hour 03 Min 29 Sec
Total Counter :	100008	Total Running Time :	3293 Hour 35 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	165	%	Manual ;	45	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		600000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	0.6	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	4.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	11.02	kN	No		
HEP :	25.0 %	10.59	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	8.65	kN	Yes 15	Times	
Ref :		8.48	kN			
LCP :	2.0 %	8.31	kN	Yes 15	Times	
LEP :	25.0 %	6.36	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	5.93	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	1.95	kN	Min :	1.85	kN	Punch NO.;	29
SD :	2.40	%	Max :	2.06	kN	Punch NO.;	49
Mean Main Pressure :	8.46	kN	Min :	8.05	kN	Punch NO.;	16
SD :	2.50	%	Max :	9.12	kN	Punch NO.;	39
Filling Depth :	Increase ;	26	Times	Decrease ;	49	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	25	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	6.90	mm	Pre Thickness ;	3.90	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	6.60	mm	Main Thickness ;	3.60	mm
Filling Depth :	8.64	mm	Current Cam :	Cam C	
Main Air pressure :	0.752	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	40	Min			
Lower Punch Remain Time (S2):	1	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0
Lower Punch :	Average ;	0.15	kN		
	Max ;	0.31	kN	Max Punch ;	20
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	100008	Tablets			
A.W.C Counter :	93744	Tablets			
HEP :	166 (0.2%)	Tablets	LEP :	36 (0.0%)	Tablets
GOOD :	93542 (99.8%)	Tablets			