

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | GEMFI<br>600MG,STKF2A00 |
| User ID :        | sj97904       | Batch NO. :          | KFM2500260              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 10 Min 02 Sec   |
| Total Counter :  | 300058        | Total Running Time : | 3295 Hour 41 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 166               | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min : | 6.50 mm              |
| Target Quantity :                  |        | 600000            | Tablets |                    | Max : | 12.00 mm             |
| Air Pressure Low Limit :           |        | 0.000             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 50                | Min     | Supply Time ;      | 0.6   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 30                | Min     | Supply Time ;      | 4.0   | Sec                  |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 2.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 1.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 10.21             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 9.82              | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 8.01              | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 7.86              | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 7.70              | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 5.89              | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 5.50              | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.54       | kN  | Min : | 1.21       | kN  | Punch NO.; | 4  |
| SD :                 | 7.80       | %   | Max : | 1.82       | kN  | Punch NO.; | 22 |
| Mean Main Pressure : | 8.17       | kN  | Min : | 7.15       | kN  | Punch NO.; | 6  |
| SD :                 | 7.40       | %   | Max : | 9.86       | kN  | Punch NO.; | 45 |
| Filling Depth :      | Increase ; | 244 | Times | Decrease ; | 201 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 25       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.90    | mm  | Pre Thickness ;  | 3.91     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.55    | mm  | Main Thickness ; | 3.55     | mm           |
| Filling Depth :               | 9.87    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.748   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 26      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 17      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |   |
|------------------|-----------|------|----|-------------|---|
| Upper Punch :    | Average ; | 0.00 | kN |             |   |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0 |
| Lower Punch :    | Average ; | 0.31 | kN |             |   |
|                  | Max ;     | 0.82 | kN | Max Punch ; | 2 |
| Ejection Force : | Average ; | 0.00 | kN |             |   |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 300058            | Tablets |       |            |         |
| A.W.C Counter : | 291971            | Tablets |       |            |         |
| HEP :           | 341 (0.1%)        | Tablets | LEP : | 193 (0.1%) | Tablets |
| GOOD :          | 291437<br>(99.8%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | GEMFI<br>600MG,STKF2A00 |
| User ID :        | sj97904       | Batch NO. :          | KFM2500260              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 10 Min 02 Sec   |
| Total Counter :  | 300058        | Total Running Time : | 3295 Hour 41 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 152               | %       | Manual ;           | 45           | RPM       |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |        | 600000            | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |        | 0.000             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 50                | Min     | Supply Time ;      | 0.6          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 30                | Min     | Supply Time ;      | 4.0          | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 2.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 1.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 11.31             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 10.87             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 8.87              | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 8.70              | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 8.52              | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 6.52              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 6.09              | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.07       | kN | Min : | 1.94       | kN | Punch NO.; | 4  |
| SD :                 | 3.50       | %  | Max : | 2.23       | kN | Punch NO.; | 22 |
| Mean Main Pressure : | 8.97       | kN | Min : | 8.41       | kN | Punch NO.; | 42 |
| SD :                 | 3.30       | %  | Max : | 9.75       | kN | Punch NO.; | 39 |
| Filling Depth :      | Increase ; | 52 | Times | Decrease ; | 82 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 23       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.90    | mm  | Pre Thickness ;  | 3.91     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.60    | mm  | Main Thickness ; | 3.61     | mm           |
| Filling Depth :               | 9.91    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.748   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 26      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 17      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |   |
|------------------|-----------|------|----|-------------|---|
| Upper Punch :    | Average ; | 0.00 | kN |             |   |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0 |
| Lower Punch :    | Average ; | 0.31 | kN |             |   |
|                  | Max ;     | 0.82 | kN | Max Punch ; | 2 |
| Ejection Force : | Average ; | 0.00 | kN |             |   |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 300058            | Tablets |       |            |         |
| A.W.C Counter : | 291971            | Tablets |       |            |         |
| HEP :           | 341 (0.1%)        | Tablets | LEP : | 193 (0.1%) | Tablets |
| GOOD :          | 291437<br>(99.8%) | Tablets |       |            |         |